

anasept[®]

ANTIMICROBIAL SKIN & WOUND CLEANSER

Ordering Information

Anasept[®] Antimicrobial Skin & Wound Cleanser

CATALOG NO.		NDC NUMBER	SIZE	CASE QUANTITY
4004C	(Dispensing Cap	67180-400-04	4 oz	24
4008C	(Dispensing Cap)	67180-400-88	8 oz	12
4008SC	(Sprayer)	67180-400-88	8 oz	12
4008TC	(Trigger Sprayer)	67180-408-88	8 oz	12
4012SC	(Trigger Sprayer)	67180-400-12	12 oz	12
4016C	(Dispensing Cap)	67180-400-16	15 oz	12


Anasept[®] Antimicrobial Skin and Wound Gel

5003G	(Tube)	67180-500-03	3 oz	12
-------	--------	--------------	------	----

Anasept[®] is a registered trademark of Anacapa Technologies. Anasept products are manufactured in the USA.



- NO KNOWN MICROBIAL RESISTANCE
- RAPID KILL
- SPORICIDAL
- SAFE AND TISSUE COMPATIBLE
- AIDS IN THE DEBRIDEMENT OF SLOUGH AND DEBRIS
- OUTSTANDING ODOR CONTROL
- HELPS REDUCE THE RISK OF INFECTION



Scanning Electron Micrograph: Staphylococcus epidermidis enlarged 27,000x



301 E. Arrow Hwy, Ste. 106
San Dimas CA 91773

Toll-Free: 800-489-2591
Direct: 909-394-7795
Fax: 909-394-9895

e-mail: anacapa@anacapa-tech.net
Website: www.anacapa-tech.net

Anasept[®] is also available as Anasept Antimicrobial Skin and Wound Gel, a clear, thick, isotonic hydrogel with long-lasting, broad-spectrum antimicrobial properties.

Anasept[®] Gel's Medicare reimbursement code is HCPCS:A6248



ANASEPT®
The clear, safe and
MOST EFFECTIVE CHOICE FOR
skin and wound antisepsis®

PRODUCT DESCRIPTION AND USES:

Anasept® is an extremely safe and gentle skin and wound cleanser with exceptionally rapid broad-spectrum bactericidal, fungicidal and virucidal properties through the action of antimicrobial sodium hypochlorite. Anasept helps in the mechanical removal of the debris and foreign material from the wound or application site. Anasept is a very pure, completely colorless, isotonic, tissue compatible solution. Anasept is stable for two-years and is free of necrotizing chemicals such as sodium hydroxide.

RAPID ACTION:

Anasept demonstrates exceptionally rapid microbicidal action. Most pathogenic organisms are killed within 2 minutes or less following application. There is no known microbial resistance to Anasept.

CLINICALLY TESTED:

Anasept is clinically proven to reduce wound bioburden levels and improve the rate of healing.*

SAFETY:

Anasept has been subjected to rigorous safety testing at an independent laboratory and shown to meet the criteria for safe use.

- Modified Primary Skin Irritation (FHSA method – 7 day exposure with repeated insult to intact and abraded skin)
- Cytotoxicity (ISO Agarose Overlay method)
- Systemic toxicity (ISO Acute Systemic Toxicity)
- ISO Sensitization Study

ENVIRONMENTALLY FRIENDLY:

Anasept does not leave any toxic residues or by-products. Anasept chemically breaks down into salt and water and is completely safe for disposal in the public sewer system.

WARNINGS:

For **External Use Only**. Discontinue use if redness or irritation develops. **Do NOT USE** in the eyes.



Scanning ElectronMicrograph: E. coli enlarged 21,000x

GENERAL DIRECTIONS FOR USE

Skin Cleansing:

- 1) Spray intended area or saturate sterile gauze with Anasept and clean skin area with a circular motion beginning at the center of the site and move outward until the selected dermal area has been thoroughly cleansed.
- 2) Air dry for 2 minutes. If preferred, allow Anasept saturated sterile gauze to remain in place as a wet dressing.

Wound Cleansing:

- 1) Debride wound, if necessary.
- 2) Spray Anasept onto entire wound bed, including the wound margin. Avoid pooling.
Alternate: Saturate sterile gauze pad with Anasept and apply to wound site.
- 3) Cover wound site with a sterile gauze or other appropriate wound dressing.
- 4) Repeat procedure once a day. Ensure that wound bed remains moist between dressing changes.

TIME KILL STUDIES

Test Organisms:

Table of Anti-microbial Activity

Pathogenic Bacteria	Initial Microorganism Count/ML	Exposure time / % Kill		
		30 seconds	1 minute	5 minutes
Escherichia coli	10 ⁸	100%	100%	100%
Staphylococcus aureus	10 ⁸	100%	100%	100%
Methicillin Resistant Staphylococcus aureus (MRSA)	10 ⁸	100%	100%	100%
Vancomycin Resistant Enterococcus faecalis (VRE)	10 ⁸	100%	100%	100%
Pseudomonas aeruginosa	10 ⁸	100%	100%	100%
Proteus mirabilis	10 ⁸	99.998%	100%	100%
Serratia marcescens	10 ⁸	100%	100%	100%
Acinetobacter baumannii	10 ⁷	–	99.96%	99.98%
Clostridium difficile	10 ⁵	100%	100%	100%
Pathogenic Fungi				
Candida albicans	10 ⁸	99.1%	99.9%	100%
Aspergillus niger	10 ⁸	99.99%	99.9999%	100%

Extremely high concentrations of pathogenic micro-organisms were exposed to Anasept over the course of precisely timed intervals in the presence of an interfering substance that simulates the organic load condition of the wound environment and is known to inhibit the action of antimicrobial agents.



TIME KILL STUDIES

Test Organisms:

Table of Sporicidal Activity

Test Substance	Initial Microorganism Count/ML	Exposure Time	Percent Reduction	Log Reduction
Clostridium difficile - spore	10 ⁶	15 minutes	99.999%	> 5.7

CATEGORIES FOR USE:

Dialysis*:
Preparation of site for Graft-Fistula Cannulation
Exit Site Dressing change for Peritoneal Dialysis
Central Line Site Preparation.

* Detailed site preparation procedures are available upon request. *Compatible with catheters used in dialysis procedures. Catheter compatibility reports available upon request.*

Topical and Wound Care:
Application to skin or wound to establish antisepsis at the site.



- **RAPID KILL**
- **NO KNOWN MICROBIAL RESISTANCE**
- **SPORICIDAL**
- **SAFE AND TISSUE COMPATIBLE**
- **AIDS IN THE DEBRIDEMENT OF NECROTIC SLOUGH AND DEBRIS**
- **OUTSTANDING ODOR CONTROL**
- **HELPS REDUCE THE RISK OF INFECTION**
- **LATEX FREE**

*J. Lindfors, A Comparison of an Antimicrobial Wound Cleanser to Normal Saline in Reduction of Biobuden and Its Effect on Wound Healing. Ostomy/Wound Management. 2004; 50 (8): 28-41.

safety

Anasept® Antimicrobial Skin and Wound Gel has been subjected to rigorous safety and toxicological evaluations to comply with FDA regulations at an independent FDA registered testing facility and shown to meet all criteria for safe use.

- Modified Skin Irritation Study (FSHA method 7 day exposure with repeated insult to intact and abraded skin)
- Cytotoxicity (USP method)
- Systemic Toxicity (USP method)
- ISO Sensitization Study.

Test reports available upon written request.

clinically tested:

Anasept® Antimicrobial Skin & Wound Cleanser, the liquid version of Anasept Antimicrobial Skin & Wound Gel is clinically proven to reduce wound bioburden levels and improve the rate of healing.*



Anasept® is also available as Anasept Antimicrobial Skin & Wound Cleanser in a wide-variety of dispensers specifically designed for skin & wound cleansing applications. Anasept Antimicrobial Skin & Wound Cleanser has all the same powerful broad spectrum antimicrobial and safety features inherent in Anasept Antimicrobial Skin & Wound Gel.

Anasept® is a registered trademark of Anacapa Technologies, Inc.

Anasept products are manufactured in the USA.

**EXTERNAL USE ONLY.
NOT FOR OPHTHALMIC USE.**

latex FREE



DIRECTIONS FOR USE:

Wound Care:

- 1) Debride wound, if necessary or cleanse wound with a wound cleanser such as Anasept® Antimicrobial Skin and Wound Cleanser.
- 2) Apply a generous amount (1/4 " to 1/2 "thick) of Anasept Antimicrobial Skin and Wound Gel to entire wound bed, including areas of undermining.
- 3) Apply a thin coating to peri-wound skin area and allow to dry.
- 4) Cover with appropriate wound dressing or covering (avoid silver and other wound dressings containing heavy metals).
- 5) Change dressing once a day. Maintain a moist wound environment between dressing changes.

NOTE: Anasept products contain sodium chloride which is not compatible with wound care products that contain silver. Silver in the presence of sodium chloride will be converted to insoluble silver chloride and become inactive.

Indwelling Vascular Catheters:

- 1) Apply sufficient quantity of Anasept Antimicrobial Skin & Wound Gel to completely cover skin area around the indwelling vascular catheter.
- 2) Cover with appropriate site dressing.

Ostomy:

- 1) Apply a thin coating of Anasept Antimicrobial Skin & Wound Gel to peri-stomal area.
- 2) Allow to dry.
- 3) Apply Ostomy appliance.

Skin Care:

- 1) Cleanse affected area with appropriate skin cleanser.
- 2) Allow to dry.
- 3) Apply a thin coating of Anasept Antimicrobial Skin & Wound Gel.
- 4) Reapply as necessary.

Ordering Information

Anasept® Antimicrobial Skin and Wound Gel

Catalog No.	NDC Number	Size	Case Quantity
5003G	67180-500-03	3oz	12
Medicare Reimbursement Code: HCPCS #A6248			

301 E. Arrow Hwy, Ste. 106

San Dimas, California 91773

Toll-Free: 800-489-2591

Tel: 909-394-7795

Fax: 909-394-9895

e-mail: anacapa@anacapa-tech.net

Website: www.anacapa-tech.net

ANASEPT® ANTIMICROBIAL SKIN AND WOUND GEL

The clear, safe
and most effective choice
for skin and wound antisepsis™

Kills all bacteria,
fungi, viruses
and spores without
harming healthy
tissue



product description:

Anasept® Antimicrobial Skin and Wound Gel is an extremely safe topical hydrogel with exceptionally rapid broad spectrum bactericidal, including the **antibiotic** resistant strains **MRSA & VRE**, fungicidal, virucidal and sporicidal properties through the action of sodium hypochlorite. There is no known microbial resistance to Anasept Antimicrobial Skin & Wound Gel.

Anasept® Antimicrobial Skin and Wound Gel is pure, completely colorless, isotonic, **non-cytotoxic**, tissue compatible viscous hydrogel. Anasept Antimicrobial Skin & Wound Gel has a **3 year shelf-life** when stored at normal room temperature up to 25° C (77° F).

time kill studies:

Extremely high concentrations of pathogenic micro-organisms were exposed to Anasept Antimicrobial Skin and Wound Gel over the course of precisely timed intervals in the presence of an interfering substance that simulates the organic load conditions of the wound environment and is know to **inhibit** the action of antimicrobial agents. Anasept Antimicrobial Skin and Wound Gel proved 100% effective against all **pathogenic micro-organisms** tested within the first ten minutes of application except for Acinetobacter baumannii where it was shown to be 99.998% effective in the same test period (see tables).

TIME KILL STUDIES

Test Organisms: Table of Antimicrobial Activity

Pathogenic Bacteria:	Initial Organism	Exposure Time/% Kill			
	Count	1 min.	3 min.	5 min.	10 min.
Escherichia coli	10 ⁷	99.25%	99.986%	99.9995%	100%
Staphylococcus aureus	10 ⁷	100%	100%	100%	100%
Methicillin Resistant Staphylococcus aureus (MRSA)	10 ⁷	100%	100%	100%	100%
Vancomycin Resistant Enterococcus faecalis (VRE)	10 ⁷	100%	100%	100%	100%
Pseudomonas aeruginosa	10 ⁷	99.996%	100%	100%	100%
Proteus mirabilis	10 ⁷	99.888%	99.998%	99.9998%	100%
Serratia marcescens	10 ⁷	100%	100%	100%	100%
Acinetobacter baumannii	10 ⁷	99.722%	99.977%	99.996%	99.998%
Clostridium difficile	10 ⁵	100%	100%	100%	100%
Pathogenic Fungi:					
Candida albicans	10 ⁶	100%	100%	100%	100%
Aspergillus niger	10 ⁶	100%	100%	100%	100%

anasept[®]
GEL

KILLS MRSA & VRE IN 60 SECONDS

indications for use:

Anasept® Gel is intended for OTC use for management of skin abrasions, minor irritations, lacerations, cuts, exit sites and intact skin.

Professional Use: Anasept® Gel is intended to be used under the supervision of a healthcare professional in the management of wounds such as stage I-IV pressure ulcers, partial & full thickness wounds, diabetic foot & leg ulcers, post surgical wounds, first & second degree burns, grafted & donor sites.

24 hour challenge test:

Anasept Antimicrobial Skin and Wound Gel was subjected to a high concentration of **pathogenic micro-organisms** (amount known to cause infection) in the presence of an interfering substance that simulates the organic load of the wound and is known to inactivate the antimicrobial agents. The duration of antimicrobial effectiveness of Anasept Antimicrobial Skin and Wound Gel was determined in a re-challenge of the original test sample with a high concentration of freshly prepared **micro-organisms** after 24 hours of initial exposure to **pathogenic micro-organisms**.

Sustained duration of action:

Anasept Antimicrobial Skin and Wound Gel was shown to maintain microbiocidal activity even after 24 hours and repeated exposure to pathogenic micro-organisms in the simulated wound environment. The gel reduced all pathogenic test organisms by more than 99% within the first fifteen minutes of repeated exposure.

TIME KILL STUDIES

Test Organism: Table of Sporicidal Activity

Test Substance	Initial Microorganism Count/ML	Exposure Time	Percent Reduction	Log Reduction
Clostridium difficile - spore	10 ⁶	15 minutes	99.986%	>4.0



TIME KILL STUDIES - 24 HOUR CHALLENGE:

Test Organisms: Table of Antimicrobial Activity

Pathogenic Bacteria:	Initial Organism Ct. / Re-challenge Organism Ct	Exposure time after re-challenge at 24 hours / % Kill		
		5 min.	10 min.	15 min.
Escherichia coli	10 ⁷ / 10 ⁷	71.25%	96.63%	99.49%
Staphylococcus aureus	10 ⁷ / 10 ⁷	95.91%	96.45%	99.16%
Methicillin Resistant Staphylococcus aureus (MRSA)	10 ⁷ / 10 ⁷	95.69%	99.38%	99.78%
Vancomycin Resistant Enterococcus faecalis (VRE)	10 ⁷ / 10 ⁷	92.8%	96.9%	99.61%
Pseudomonas aeruginosa	10 ⁷ / 10 ⁷	84.35%	98%	99.88%
Proteus mirabilis	10 ⁷ / 10 ⁷	67.14%	97.71%	99.74%
Serratia marcescens	10 ⁷ / 10 ⁷	96%	99.36%	99.94%
Acinetobacter baumannii	10 ⁷ / 10 ⁷	13.64%	85.25%	99.25%
Pathogenic Fungi:				
Candida albicans	10 ⁶ / 10 ⁶	98.89%	99.99%	99.9996%
Mix of all above including Candida albicans	10 ⁷ / 10 ⁷	88.75%	97.31%	99.8%

*J. Lindfors. A Comparison of an Antimicrobial Wound Cleanser to Normal Saline in Reduction of Bioburden and its Effect on Wound Healing. Ostomy / Wound Management 2004; 50 (8): 28-41.

GET IN THE FRESH ZONE SANI-ZONE™



The most effective way to
get rid of those nasty
odors in the air or in the
ostomy pouch



800-489-2591
301 East Arrow Highway
Suite 106
San Dimas, CA 91773

*Safe - certified non-toxic by the Consumer
Product Testing Company, Inc.*

Lasts for hours! Extraordinary odor eliminating capability. Leaves a squeaky clean background smell. The Best product you have ever used!

**Sani-Zone™ Maximum
Strength Odor
Eliminator Spray**



1008A

1002A

One or two sprays in an average size room eliminates odors and leaves a squeaky clean fragrance

*Outstanding odor control!
Unmatched in effectiveness!*

**Sani-Zone™ Ostomy
Appliance Deodorant**



1008-OD

1002-OD

Just five drops in an ostomy pouch stops odors for hours.

Powerful and fast-acting, it stands up against the most embarrassing odors.

*Medicare Reimbursement HCPCS Code
for Sani-Zone Ostomy Appliance
Deodorant: A4394*

GET IN THE FRESH ZONE SANI-ZONE™



The most effective way to
get rid of those nasty
odors in the air or in the
ostomy pouch



800-489-2591
301 East Arrow Highway
Suite 106
San Dimas, CA 91773

*Safe - certified non-toxic by the Consumer
Product Testing Company, Inc.*

Lasts for hours! Extraordinary odor eliminating capability. Leaves a squeaky clean background smell. The Best product you have ever used!

**Sani-Zone™ Maximum
Strength Odor
Eliminator Spray**



1008A

1002A

One or two sprays in an average size room eliminates odors and leaves a squeaky clean fragrance

*Outstanding odor control!
Unmatched in effectiveness!*

**Sani-Zone™ Ostomy
Appliance Deodorant**



1008-OD

1002-OD

Just five drops in an ostomy pouch stops odors for hours.

Powerful and fast-acting, it stands up against the most embarrassing odors.

*Medicare Reimbursement HCPCS Code
for Sani-Zone Ostomy Appliance
Deodorant: A4394*

Product Description:

Silver-Sept® Silver Antimicrobial Skin & Wound Gel is a clear, amorphous hydrogel wound dressing that helps to maintain a moist wound environment that is conducive to healing. Silver-Sept will absorb a mild amount of wound exudate. Silver-Sept will not stain or discolor tissue.

Silver-Sept functions as a long-lasting, antimicrobial barrier by inhibiting the growth of bacteria including the antibiotic re-sistant strains: MRSA & VRE, as well as, fungi such as: Candida albicans and Aspergillus niger.

Biocompatibility and Safety:

Silver-Sept has been subjected to rigorous biocompatibility and safety testing at independent testing laboratories and shown to be non-irritating, non-sensitizing and non-cytotoxic and meets the highest standards for safe use. Testing included:

- Cytotoxicity Studies
- ISO Skin Irritation Studies
- ISO Sensitization Studies
- USP and ISO Modified Systemic Toxicity Studies

Sustained Duration of Action:

Silver-Sept has been shown in antimicrobial effectiveness studies (see table 2) to maintain its high level of antimicrobial barrier properties for up to 3 days.



Antimicrobial Activity

TIME KILL STUDIES: The time kill studies for Silver-Sept® Silver Antimicrobial Skin and Wound Gel, were conducted in the presence of an interfering substance that simulates the organic load conditions of the wound environment and is known to inhibit the action of antimicrobial agents.

Test Organisms:		Table 1 – Antimicrobial Activity			
Bacteria:	INITIAL ORGANISM COUNT	EXPOSURE TIME / % KILL			
		1 hour	2 hours	1 day	
Escherichia coli	10 ⁷	99.99985%	100 %	100 %	
Staphylococcus aureus	10 ⁷	99.77%	100 %	100 %	
Methicillin Resistant Staphylococcus aureus (MRSA)	10 ⁷	98.32%	100 %	100 %	
Vancomycin Resistant Enterococcus faecalis (VRE)	10 ⁷	98.27%	100 %	100 %	
Pseudomonas aeruginosa	10 ⁷	99.9996%	100 %	100 %	
Proteus mirabilis	10 ⁷	99.9998 %	100 %	100 %	
Serratia marcescens	10 ⁷	99.9538%	100 %	100 %	
Fungi:					
Candida albicans	10 ⁶	99.1 %	99.9 %	100 %	
Aspergillus niger	10 ⁶	99.99 %	99.9999 %	100 %	

TIME KILL STUDIES: The time kill studies for Silver-Sept® Silver Antimicrobial Skin and Wound Gel, were conducted by inoculating high level of test organisms, repeatedly for three consecutive days.

Test Organisms:		Table 2 – Repeated Exposure Time Kill Study					
Bacteria:		EXPOSURE TIME / % KILL					
		2nd day			4 days after 3rd		
		1 day after	2nd	1 day after	3rd	3rd day, 1 day	
		initial test	inoculation	2nd	inoculation	after 3rd	
		inoculum	inoculation after day 1	inoculation	after day 2	inoculation	
Escherichia coli	10 ⁷	100 %	10 ⁷	100 %	10 ⁷	100 %	100 %
Staphylococcus aureus	10 ⁷	100 %	10 ⁷	100 %	10 ⁷	100 %	100 %
Methicillin Resistant Staphylococcus aureus (MRSA)	10 ⁷	100 %	10 ⁷	100 %	10 ⁷	100 %	100 %
Pseudomonas aeruginosa	10 ⁷	100 %	10 ⁷	100 %	10 ⁷	100 %	100 %
Serratia marcescens	10 ⁷	100 %	10 ⁷	100 %	10 ⁷	100 %	100%
Fungi:							
Candida albicans	10 ⁶	100 %	10 ⁶	100 %	10 ⁶	100 %	100%
Aspergillus niger	10 ⁶	98.667 %	10 ⁶	98.333%	10 ⁶	99.5%	100%

Indications:

Silver-Sept® is intended for OTC use for Abrasions and Lacerations and under the supervision of a healthcare professional in the management of:

- Stage I- IV Pressure Ulcers
- Partial and Full Thickness wounds
- Diabetic Foot and Leg Ulcers
- 1st & 2nd degree burns
- Grafted and Donor Sites

Directions for Use:

Wound Care:

- Cleanse or debride wound as necessary.
- Apply a generous amount of Silver-Sept Gel directly into wound bed (1/8” to 3/16”).
- Cover with a sterile gauze or other appropriate secondary dressing and secure in place.
- Maintain a moist wound environment between dressing changes.

Note: Silver-Sept Silver Antimicrobial Skin & Wound Gel may remain in the wound bed for up to 3 days. More frequent dressing changes may be required dependent upon the amount of wound exudate present and the condition of the secondary dressing.

Ostomy:

- Apply a thin coating of Silver-Sept to peristomal area.
- Allow to dry.
- Apply Ostomy appliance.

Skin Care:

- Cleanse affected area with appropriate skin cleanser.
- Allow to dry.
- Apply a thin coating of Silver-Sept.
- Repeat as necessary.

EXTERNAL USE ONLY. NOT FOR OPHTHALMIC USE
latexFREE

Silver-Sept® is a trademark of Anacapa Technologies, Inc.
US & International Patents Pending

Ordering Information

Silver-Sept® Silver Antimicrobial Skin & Wound Gel ordering information:

Catalog no.	Description	Packaging
3015 S	1.5 oz tube	12/cs
3003 S	3 oz. tube	12/cs



Medicare Reimbursement Code: HCPCS A6248



Silver-Sept® is manufactured in the USA by:
Anacapa Technologies, Inc.
San Dimas, CA
Tel: (800) 489-2591
Fax: (909) 394-9895
e-mail: anacapa@anacapa-tech.net
Website: www.anacapa-tech.net

