

# Bilirubin Confirmation Made Easy!



## Positive Results = Blue or Purple Color<sup>1</sup>


The test is positive if a blue or purple color develops on the test mat at 60 seconds. This indicates that bilirubin is present. A slight pink or red color should be ignored.

## Negative Results = No Blue or Purple Color

The test is negative if no blue or purple color develops on the test mat within 60 seconds.

## AimTab™ Bilirubin Tablets

Same test procedure as Ictotest®, a registered trademark of Siemens Healthcare.



CLIA  
WAIVED

- AimTab™ Bilirubin Tablets, 100/bottle, REF #13200
- Provides quick & easy urine bilirubin confirmatory results
- Detects 0.05 to 1.0 mg/dL (0.8 to 1.7  $\mu$ mol/L)
- Confirms "urine strip" results
- Eliminates the need to perform blood test for bilirubin
- CAP / CLSI formatted procedures available via fax or email - just add your name!
- Store at room temperature.

Tests you Want,  
Accuracy you Need!

<sup>1</sup> See Package Insert for complete procedure.

©2014 Germaine Laboratories, Inc., CSM 0800-19, Rev. 2-14

Distributed by:

**GERMAINE**<sup>®</sup>  
LABORATORIES, INC.

[www.germainelabs.com](http://www.germainelabs.com)

800.854.8446

210.692.4192

210.692.4198 Fax

