

Your Intelligent, Total Fluid Delivery System

Syringe and Disposable Kits

#### **Designed for Optimal Workflow**

 Quick, Alignment-Free Attachment and Detachment

Saves time and maintains aseptic conditions

Automatic Plunger Sensing System
 Detects and reacts to syringe state

Footless Plunger
 Ensures correct piston-plunger alignment every time

# Medrad Disposable Features from Bayer

- Ultra-clear barrel permits visualization of air bubbles
- Patented Medrad® FluiDot® Indicators allow operators to detect the presence of fluid
- Prime tube captures fluid overflow to minimize drip & mess during priming

### **Additional Options**

- Sterile spike minimizes handling of contrast bottle
- Transfer sets economically enhance fluid loading efficiency





**MEDRaD®** Stellant®


CT Injection System



Swahhahle-valve

transfer set

50 per box



## How to Order Quality Medrad® Stellant® Disposable Products from Bayer

Our syringe kit selection allows you to choose the optimum configurations for your procedures. Talk to your Bayer Representative for more information.

#### **New! Order disposables online**

The myORDERS website makes it easy to order Bayer disposables and keep track of your orders online. We take care of your orders so you can take care of your patients.





Syringe + 60" low-pressure

tubing with prime tube +

quick-fill tube

50 per box

Pharmaceuticals Division

Bayer HealthCare LLC 100 Bayer Boulevard P.O. Box 915 Whippany, NJ 07981 www.radiologysolutions.bayer.com

Bayer Medical Care Inc. 1 Bayer Drive Indianola, PA 15051 USA Customer Service/Orders 1-800-633-7231 Customer Service Fax 1-412-767-4120 Bayer reserves the right to modify the specifications and features described herein, or discontinue manufacture of the product described at any time without prior notice or obligation. Please contact your authorized Bayer representative for the most current information.

Bayer, the Bayer Cross, Medrad®, Stellant®, and FluiDot® are trademarks of the Bayer group of companies.

© 2011. 2014. 2016 Bayer

PP-STL-D-US-0028 November 2016