

United States

Specialty Products Catalog

The Fresenius Kabi line of specialty products is designed to help operators process blood with accuracy and efficiency.





Directory

Fresenius Kabi USA, LLC Three Corporate Drive Lake Zurich, Illinois 60047

Main

800-333-6925

- · Customer Service
- Technical Services
- Quality Assurance
- Product Information
- · Accounts Receivable
- · Accounts Payable

For more information, please contact us at productsupport@fresenius-kabi.com or 800-333-6925.

Refer to the product operator's manual, instructions for use, and other associated labeling and information for a complete list of warnings, cautions, contraindications, and other pertinent information associated with the use of the products in this catalog.

Table of Contents

Specialty Products	Page(s)
Processing Disposables	1
Transfusion/Testing/Therapeutic	1-2
CompoMat G5 & CompoMaster NET System	3
Transfer-Pack Containers	4
Storage/Freezing	5-6
Leukoreduction Filters	7-8
Sealers	9
CompoDock DM & DockMaster NET System	9-10
Processing Hardware	10
CompoLab TS	11
Motorized Blood Donor Chair	11
CompoCool	12
CompoGuard & DonationMaster NET System	13-14
Platelet Bacterial Testing	15
Index	16

Processing Disposables

Code	ltem	Quantity
4C2405	Fenwal Sampling Site Coupler Natural rubber sleeve stopper allows sampling with a 22 gauge or smaller needle.	96 units/case
4R4350 OTERUMO' SCO WAFERS SACRIFUS BOOK OF THE SA	Terumo Welding Wafers - Disposable Needed to seal, cut and connect tube to another tube while keeping internal sterility.	140 units/case
4R4409	Plasma Container Carton Provides added protection during storage.	144 units/case
4R4423H PAT PEND. MADE IN UBA	Side Clamps (bulk pack) Provides a temporary closure to the tubing.	500 units/case
9024531	Scissor Clamp Hemostat, used to clamp tubing off.	1 each

Transfusion/Testing/Therapeutic

Code	Item	Quantity
A100645	Samplok Safety Vacuum Tube Holder Multi-Sample tube holder used to collect blood samples into vacuum test tubes and long neck blood culture bottles.	425 units/case
A100750	Samplok Sampling Kit Sampling and transfer kit for bacterial testing; 10 mL single unit chamber.	110 units/case
4C2160	Fenwal Blood Component Recipient Set with Standard Blood Filter and Female Luer Adapter Standard 170-260 micron filter, transfuses all components.	48 units/case

Transfusion/Testing/Therapeutic

Code	Item	Quantity
4C2161	Fenwal Y-Type Blood Component Recipient Set with Standard Blood Filter and Luer Adapter Standard 170-260 micron filter with Y-set, transfuses all components.	48 units/case
4C2223	Blood Component Infusion Set Blood component infusion set with female luer connector for syringe, transfuses all components, distal needle/catheter adapter contains 170-260 micron blood component filter, 13" in length.	48 units/case
4C2240	2 Units Fenwal Plasma Transfer Set Plasma transfer set with coupler and male luer lock, 30" long.	96 units/case
4C2243	2 Units Fenwal Plasma Transfer Set Plasma transfer set with two couplers, 30" long.	96 units/case
4C2498H	Y-Type Connector Set with Luer Lock Commonly used for autologous transfusions with Fenwal BPU 4RO499.	48 units/case
4R5125	Needle Guard Post collection needle stick protection device.	1,000 units/case
4R5126	Hematype Segment Device Single use sampling device designed to help safely access blood from a segment.	2,000 units/case
4R5128	Hematype Segment Device Single use sampling device designed to help safely access blood from a segment.	1,000 units/case

CompoMat G5 & CompoMaster NET System

Automated whole blood and blood component separator (post-centrifugation). Includes 3 independent presses, up to 6 independent clamps/sealers, 2 detector types with up to 13 detector locations, 3 scales, 2 suspension pins; internal memory of 1000-2000 preparations.

9025561 CompoMat G5



Code	Item	Quantity
9025521	PLS Weighing Scale Weighing scale designed to provide tared or untared weight of a plasma blood product in BPU systems. Includes one scale, one clamp/sealer, one deaeration plate, one detector. One PLS Weighing Scale suggested per CompoMat G5.	1 unit
9025511	RBC Weighing Scale Weighing scale to provide tared or untared weights. One RBC Weighing Scale recommended per CompoMat G5 for processes using top & bottom BPU systems to weigh red cell blood products; includes a passive metal insert to assist in manual cannula breaking.	1 unit
9025571	CompoMaster NET Data management software dedicated to the CompoMat G5. Designed to set up programs on and capture data after each process from the CompoMat G5. Software is provided on a CD-Rom; computer hardware not included; one CompoMaster NET required per computer set up; one CompoMaster NET can manage multiple CompoMat G5 devices.	1 unit
9025611	Scanner Barcode scanner dedicated to CompoMat G5 for data capture during whole blood or blood component processing, such as operator barcode, donation ID number, pass/fail result codes, centrifuge number, etc. One scanner suggested per CompoMat G5.	1 unit
M684151	WLAN Access Point Dedicated electrical device designed to send and receive data wirelessly between a computer and the CompoMat G5.	1 unit
M689401	WLAN USB Stick USB stick designed to send data wirelessly from a single CompoMat G5 device to a WLAN Access Point and connected computer.	1 unit
M684201	LAN Hub 5-Port Dedicated electrical device designed to send and receive wired data from one computer to multiple CompoMat G5 devices. Up to 5 connection ports provided. 5m cable for CompoMat G5 not included.	1 unit
M684211	LAN Hub 24-Port Dedicated electrical device designed to send and receive wired data from one computer to multiple CompoMat G5 devices. Up to 24 connection ports provided. 5m cable for CompoMat G5 not included.	1 unit

Transfer-Pack Containers

Transfer-Pack containers in multiple configurations to address your needs. Dispensing cartons available on most Transfer-Packs.

Code	Item	Quantity
4R2001	Transfer-Pack container with Coupler – 150 mL 8 packs of 12 in each dispensing carton.	96 units/case
4R2004	Quadruple Transfer-Pack container with Coupler – 150 mL 8 packs of 6 in each box (not available in dispensing carton).	48 units/case
4R2014	Transfer-Pack container with Coupler – 300 mL 8 packs of 12 in each dispensing carton.	96 units/case
4R2023	Transfer-Pack container with Coupler – 600 mL 8 packs of 12 in each dispensing carton.	96 units/case
4R2024	Transfer-Pack container with Male Luer – 600 mL 8 packs of 12 in each dispensing carton. Empty bag with male luer lock.	96 units/case
4R2027	Transfer-Pack container with 8 Couplers – 600 mL Eight couplers allow for pooling of products from multiple units.	48 units/case
4R2032	Transfer-Pack container with Coupler – 1000 mL 8 packs of 12 in each dispensing carton.	96 units/case
4R2255P	Quadruple Transfer-Pack container with Male Luer - 1000 mL Intended for plasma collection with male luer adapter, one 1000 mL and three 300 mL containers (not available in dispensing carton).	30 units/case


Storage/Freezing

Code		Item	Quantity
4A7831	Fig. 1 Sec. 1 Se	Glycerolyte 57 Buffered Glycerolyte - USP 400 mL bottle Buffered Glycerine solution for freezing of red blood cells.	12 units/case
4A7833	The state of the s	Glycerolyte 57 Buffered Glycerolyte - USP 500 mL bottle Buffered Glycerine solution for freezing of red blood cells.	12 units/case
4B7870X	Transfer of the second	1.6% Sodium Chloride Solution - USP 1000 mL Wash solution for glycerolyte.	14 units/case
4B7871	TOTAL CONTROL OF THE PARTY OF T	1.6% Sodium Chloride Solution – USP 2000 mL Wash solution for glycerolyte, larger container.	6 units/case
4B7874Q	CARCO MILES	12% Sodium Chloride – USP 150 mL Wash solution for glycerolyte.	36 units/case
4B7877X	Table 1 Control of the Control of th	0.2% Dextrose and 0.9% Sodium Chloride Solution – USP 1000 mL Wash solution for glycerolyte.	14 units/case
4B7878	The second secon	0.2% Dextrose and 0.9% Sodium Chloride Solution – USP 2000 mL Wash solution for glycerolyte, larger container.	6 units/case

Storage/Freezing

Code	Item	Quantity
4R2422	Red Blood Cell Freezing Pack with Drip Chamber and Coupler - 2000 mL Commonly used to glycerolize blood, large size allows for the RBC's glycerolyte and wash solution.	48 units/case
4R4413P	Blood-Pack Unit Holder – 2000 mL Provides added protection to Blood-Pack unit during storage.	150 units/case
4R4420	Cryoprecipitate Carton Provides added protection to cryoprecipitate during storage.	144 units/case
4R4424	Plastic Centrifuge Bag Holds BPU together during centrifugation.	1,000 units/case


Leukoreduction Filters



for Lab Sets

Code	Item	Quantity
4C2300	Sepacell R-500II Pre-Storage Leukocyte Reduction Set for Red Cells Sepacell R-500II pre-storage leukocyte reduction set with spike for use with one unit of red cells or whole blood in additive solution.	18 units/case
4C2468	Sepacell R-500II Leukocyte Reduction Filter for Red Blood Cells Sepacell R-500II leukocyte reduction filter with no pre-attached tubing for use with 1 unit of red cells.	18 units/case
4C4300	Leukoconnect Pre-Storage Sterile Connect Leukoreduction Set with Sepacell RS-2000 Filter Leukoconnect pre-storage sterile connect leukocyte reduction set with Sepacell RS-2000 filter for red blood cells in additive solution.	30 units/case

Leukoreduction Filters



for Bedside Administration

Code	Item	Quantity
4C4502	Sepacell PLS-10A Leukocyte Reduction Y-Administration Set Sepacell PLS-10A leukocyte reduction filter with Y-administration set below filter for use with one single donor or up to 10 pooled random donor platelet units.	18 units/case
4C2486	Sepacell R-500II Leukocyte Reduction Filter with Administration Set for Red Blood Cells Sepacell R-500II leukocyte reduction filter with administration set for use with 1 unit of red cells.	18 units/case
4C2491	Sepacell R-500II Leukocyte Reduction Filter with Y-Administration Set for Red Blood Cells Sepacell R-500II leukocyte reduction filter with Y-administration set for use with 1 unit of red cells.	18 units/case
4C2493	Sepacell PLS-5A Leukocyte Reduction Filter for Platelets Sepacell PLS-5A leukocyte reduction filter with no pre-attached tubing for use with one single donor or up to 6 pooled random donor platelet units.	18 units/case

Sealers

A line of tube sealers, including portable and bench-top units, designed for your sealing needs.

Code	Item	Quantity
9027011	CompoSeal Mobilea II Set Battery operated, radio frequency (RF) sealer for medical PVC tubing; set includes yellow handheld sealer, Powerpack, Charger and Transport case.	1 unit
M685181	Shoulder Pouch for CompoSeal Mobilea II Black nylon pouch holds CompoSeal Mobilea II Powerpack and handheld sealer; includes shoulder strap and belt loop (black).	1 unit
9025661	CompoSeal Universal Robust, table top, radio frequency (RF) sealer for medical PVC tubing up to 6mm OD; includes autotuning to adapt sealer to different tube sizes (purple).	1 unit
9024031	Handheld Sealer Handheld sealer with RF cable doubles CompoSeal Universal sealing capacity (purple).	1 unit
9028811	Transport Case for CompoSeal Universal Protective case to transport CompoSeal Universal (gray).	1 unit
9029700	Handsealer Splash Guard Can be used with Mobilea II, CompoGuard handheld sealer, and CompoSeal Universal optional handheld sealer.	6 units

CompoDock DM & DockMaster NET System

Sterile Connection Device

CompoDock DM & DockMaster NET System is a sterile connection device for medical PVC tubing; does not require wafers; includes counter (100-count) for training and validation.



Code	Item	Quantity
9028721	CompoDock DM Waferless sterile connection device for medical PVC tubing; includes counter (100-count) for training and validation.	1 unit
9028681	CompoDock Counter 2500 Disposable counter for up to 2500 docks and device QC testing.	1 unit

CompoDock DM & DockMaster NET System

Sterile Connection Device (Continued)



Code	Item	Quantity
9028691	CompoDock Counter 1000 Disposable counter for up to 1000 docks and device QC testing.	1 unit
9028701	CompoDock Counter 500 Disposable counter for up to 500 docks and device QC testing.	1 unit
9028921	DockMaster NET Data management software dedicated to the CompoDock DM; designed to set up programs on and capture data after each process from the CompoDock DM. Software is provided on a CD-Rom. Computer hardware not included. One DockMaster NET required per computer set up. One DockMaster NET can manage multiple CompoDock DM devices.	1 unit
9025611	Scanner Barcode scanner dedicated to CompoDock DM for data capture during the docking process, such as operator barcode, donation ID numbers, pass/fail result codes, leukoreduction filter or aliquot system batch numbers, etc.	1 unit

Processing Hardware

Code	Item	Quantity
4R4414	Manual Plasma Extractor Provides constant pressure from bottom up to extract plasma/ platelets from whole blood.	1 unit
4R4418	Clips, Handsealer Create a hermetic seal for the tubing (use with 4R4454).	1,000 units/case
4R4454	Handsealer Create a hermetic seal for the tubing.	1 unit
6R4452	Tube Stripper Tube stripper used to strip undiluted blood from donor tubing into collection container.	1 each
FESRD0500	Leuko Pole with Blood-Pack Unit Holder, Tray Allows for filtration of blood, lower base makes the pole more stable, includes lower hooks.	1 unit

CompoLab TS

Hemoglobin Analyzer

CompoLab TS is a hemoglobin analyzer used prior to whole blood or component donations.



Code	Item	Quantity
9033300	CompoLab TS Hemoglobin analyzer; includes transport case, power adapter and USB cable.	1 unit
9033310	CompoLab TS Cuvettes 500 cuvettes per case (provided in 5 foil pouches with 100 cuvettes each), single use disposable for CompoLab TS.	1 unit
908.0014	Control Low-Medium-High QC controls for low, medium and high hemoglogin. Each pack contains 1 small bottle each of low, medium and high.	1 unit

Motorized Blood Donor Chair



Code	Item	Quantity
4R4211	Motorized Blood Donor Chair, Decorative Design This Donor Chair is fully upholstered and designed to provide comfort and features variable positioning for either arm, as well as reclining and upright body positions. The donor chair is available in a wide variety of colors.	1 unit
5R5211	Motorized Blood Donor Chair, Smooth Design This Donor Chair is fully upholstered and designed to provide comfort and features variable positioning for either arm, as well as reclining and upright body positions. The donor chair is available in a wide variety of colors.	1 unit

CompoCool

Whole Blood Transport System

CompoCool is a transport system designed to cool whole blood directly after collection to room temperature.





9024251	and	90	2126

Code	Item	Quantity
9028751	CompoCool WB Cooling Cartridge Cartridge with temperature indicator for rapid cooling of up to 6 whole blood units to room temperature; dedicated cartridge for CompoCool WB insulating box.	1 unit
9028841	CompoCool WB Insulating Box Insulating box with lid designed to hold CompoCool WB cartridge and up to 6 units of whole blood.	1 unit
9028831	Trolley for CompoCool WB Cart to stack and transport up to 4 CompoCool WB insulation boxes with cooling cartridges and 24 whole blood units.	1 unit
M658391	Box Seals Seals designed to lock CompoCool WB insulating box allowing blood center to check box was not tampered with during transportation; 2000 seals per bag.	1 unit
9024251	CompoCool Cooling Cartridge Cartridge for rapid cooling of up to 8 whole blood units to room temperature; dedicated cartridge for CompoCool crate.	4 units
9024261	CompoCool Crate Crate designed to hold CompoCool cartridge and up to 8 units of whole blood; crates easily stack for both transport and storage.	1 unit
9024271	Trolley Base Trolley base (360 deg wheels) designed to stack and transport up to 4 CompoCool crates, cartridges and 32 whole blood units.	1 unit
9024541	Trolley Push Bar Push bar for trolley base for transportation.	1 unit

CompoGuard & DonationMaster NET System

Whole Blood Collection System

CompoGuard & DonationMaster NET System is a mixing scale and data management software system intended for the collection of whole blood at fixed sites and in mobile donor operations.



Code	Item	Quantity
9029661	CompoGuard COMPLETE Includes CompoGuard mixing scale, adapter, battery, tray, filter holders, gooseneck, donor display, optenna, scanner holder, handsealer holder, sealing generator, and DECT* option.	1 unit
9029651	CompoGuard DATA Includes CompoGuard mixing scale, adapter, tray, filter holders, and DECT* option (battery sold separately).	1 unit
9029641	CompoGuard BASIC Includes CompoGuard mixing scale, tray, adapter, filter holders (battery sold separately).	1 unit
9029331	DonationMaster NET Data management software dedicated to the CompoGuard; designed to set up programs on and capture data after each whole blood collection on the CompoGuard; software provided on CD-Rom; computer hardware not included; one DonationMaster NET required per computer set up; one DonationMaster NET can manage multiple CompoGuard devices.	1 unit
9029671	Barcode Scanner Barcode scanner with auto-sense function for hands-free scanning of barcodes; dedicated to CompoGuard for data capture during whole blood collections, such as operator barcode, donor and/or donation ID number, pass/fail result codes, sample tube numbers, etc.; one scanner required per CompoGuard.	1 unit
9029361	RF Converter (Fixed Part) Computer receiver for radio transmission of data signal via secured band width (DECT*) to and from CompoGuard mixing scales; one Fixed Part required per computer set up.	1 unit

^{*}DECT = Digital Enhanced Cordless Telecommunication technology

CompoGuard & DonationMaster NET System

Whole Blood Collection System (Continued)

Code	Item	Quantity
9029491	Multi-Battery Charger for CompoGuard Recharges up to 4 batteries.	1 unit
9028601	Network Converter (RS485) Computer receiver for wired data signal transmission to and from CompoGuard mixing scales; one network converter required per computer set up.	1 unit
M687051SP	Rechargeable Battery Includes battery holder.	1 unit
9029531	Gooseneck Includes donor display and optenna.	1 unit
M689751SP	Handheld Sealer White handheld sealer/wand for CompoGuard Complete with RF cable. (included with CompoGuard Complete)	1 unit
9029481	Transport Case Includes an integrated charging adapter to recharge the battery without unpacking the device.	1 unit
M636702	Y-Piece For network wired data transfer between DonationMaster NET and multiple CompoGuard devices (one y-piece per CompoGuard device).	1 unit
M636672	Network Cable For network wired data transfer between DonationMaster NET and a CompoGuard device; when using y-piece to chain multiple devices to each other, one network cable is required per CompoGuard.	1 unit
M636301	Modification Kit – Sealer Upgrade** Includes handsealer, generator and RF cable; designed for upgrade of CompoGuard BASIC or CompoGuard DATA.	1 unit
M690041	Modification Kit – DECT* Upgrade** Includes protected radio band (DECT) for wireless data transmission; designed for upgrade of CompoGuard BASIC.	1 unit
9028391	Software License Software license key required when upgrading a CompoGuard BASIC to include data management (DATA or COMPLETE). One software license activates the data management option for one CompoGuard.	1 unit
9029700	Handsealer Splash Guard Can be used with Mobilea II, CompoGuard handheld sealer, and CompoSeal Universal optional handheld sealer.	6 units

^{*}DECT = Digital Enhanced Cordless Telecommunication technology

^{**}Upgrades require Depot Service and related charges.

Platelet Bacterial Testing



Code	ltem	Quantity
P100	Platelet PGD Test (without Controls) (1) Box of 100 Platelet PGD Test Devices (P100T)	100 tests
	(1) Box of Reagents (6 Bottles, 2 each of Reagents, 1, 2, & 3) for 100 Tests (P100R)	
	(1) Box of Disposable Pipettes & Microfuge Tubes for 100 Tests (P100D)	
P20	Platelet PGD Test (without Controls) (1) Box of 20 Platelet PGD Test Devices & Disposable Pipettes & Microfuge Tubes (P20T)	20 tests
	(1) Box of Reagents (3 Bottles, 1 each of Reagents 1, 2 & 3) for 20 Tests (P2OR)	
P30C	Platelet PGD Controls (1) Box of Controls (2 Bottles, 1 Positive & 1 Negative)	1 set
P00614	PGD Sample Processing Tray (1) Each Sample Processing Tray	1 each

Index

Code	Page	Code	Page	Code	Page	Code	Page
A100645	1	4C2240	2	4R4454	10	9028841	12
A100750	1	4C2243	2	4R5125	2	9028921	10
FESRD0500	10	4C23OO	7	4R5126	2	9029331	13
M636301	14	4C24O5	1	4R5128	2	9029361	13
M636672	14	4C2468	7	5R5211	11	9029481	14
M636702	14	4C2486	8	6R4452	10	9029491	14
M658391	12	4C2491	8	90B.0014	11	9029531	14
M684151	3	4C2493	8	9024031	9	9029641	13
M684201	3	4C2498H	2	9024251	12	9029651	13
M684211	3	4C4300	7	9024261	12	9029661	13
M685181	9	4C4502	8	9024271	12	9029671	13
M687051SP	14	4R2001	4	9024531	1	9029700	9, 14
M689401	3	4R2004	4	9024541	12	9033300	11
M689751SP	14	4R2014	4	9025511	3	9033310	11
M690041	14	4R2023	4	9025521	3		
P00614	15	4R2024	4	9025561	3		
P100	15	4R2027	4	9025571	3		
P20	15	4R2O32	4	9025611	3, 10		
P30C	15	4R2255P	4	9025661	9		
4A7831	5	4R2422	6	9027011	9		
4A7833	5	4R4211	11	9028391	14		
4B7870X	5	4R4350	1	9028601	14		
4B7871	5	4R4409	1	9028681	9		
4B7874Q	5	4R4413P	6	9028691	10		
4B7877X	5	4R4414	10	9028701	10		
4B7878	5	4R4418	10	9028721	9		
4C2160	1	4R4420	6	9028751	12		
4C2161	2	4R4423H	1	9028811	9		
4C2223	2	4R4424	6	9028831	12		

CompoDock and CompoMat are trademarks of Fresenius Kabi.

PGD is a trademark of Verax Biomedical Incorporated.

Sepacell is a trademark of Asahi Kasei Medical Co.

SampLok is a trademark of Noble House Group Pty Ltd.

Over time this information may be subject to change based upon modifications in product manufacturing or changes in blood banking or clinical practice. Further inquiries may be directed to your local representative or our hotline at 1-800-937-5060



