

RECEPTAL® Floo	orstands and Brackets. Most are stainless steel - Latex-Free	
List No.	Description	Case Size
43425-05	Short canister support for wall bracket mounting	5
List No.	Description	Case Size
43426-05	Long canister support for wall bracket mounting	5
List No.	Description	Case Size
43428-05	Long canister support with mounting plate for vacuum gauge	1
List No.	Description	Case Size
43429-01	Long tandem canister support with single slide for wall or IV bracket mounting	1
List No.	Description	Case Size
43431-01	Wire floor stand	1
List No.	Description	Case Size
43432-01	Tandem floor stand with conductive casters	1
List No.	Description	Case Size
43447-01	In-line quad floor stand with conductive casters	1

RECEPTAL® Acc	essories – Latex-Free	
List No.	Description	Case Size
43316-01	Double 1/4" male Receptal® taper (for patient port to tubing connection)	100
List No.	Description	Case Size
43448-25	1/4" male x 1/4" female Receptal® taper	25
List No.	Description	Case Size
43453-25	3/8" male x 1/4" female Receptal® taper	25
List No.	Description	Case Size
43473-25	Double 1/4" male Receptal® taper	25
List No.	Description	Case Size
43494-25	1/4" male x 1/4" female Receptal® elbow	25

List No. 44300-05	Description Double 1/4" male "T" with "O" ring and screw	Case Size 5
List No. 44310–25	Description 1/4" male x 3/8" female Receptal® taper	Case Size
List No.	Description	0 6:
44311-25	1/4" male x 3/8" female orange elbow	Case Size 25

EZE-VAC® Floorstands and Brackets. Most are stainless steel – Latex-Free			
List No.	Description	Case Size	
43467-12	Gomco mounting bracket for EZE-VAC® canisters	12	
List No.	Description	Case Size	
43469-12	Ring canister support for EZE-VAC® canisters	12	

EZE-VAC® Accessories - Latex-Free			
List No. 43313-25	Description 1/4" male x 3/8" female EZE-VAC® elbow	Case Size 25	
List No. 43317-25	Description EZE-VAC® 18" transfer tubing	Case Size	

RECEPTAL® —1000 Series (Liner capacity is approximate) – Latex–Free			
List No. 43449-01	Description 1000 mL Receptal® eanister (for use with 1000 mL liners)	Case Size 10	
List No. 43056-01	Description 1000 mL liner with VAC-GARD® and non-pour spout lid	Case Size 50	

RECEPTAL® —1500 Series (Liner capacity is approximate) – Latex-Free			
List No.	Description	Case Size	
43423-01	1500 mL Receptal® canister (for use with 1500 mL liners)	10	
List No.	Description	Case Size	
43025-01	1500 mL liner with pour spout lid	50	
List No.	Description	Case Size	
43044-01	1500 mL liner with VAC-GARD® and pour spout lid	50	
List No.	Description	Case Size	
43038-05	1500 mL liner with VAC-GARD® and 144" patient tubing and pour spout lid	25	
List No.	Description	Case Size	
43064-05	2000 mL liner with VAC-GARD® and anti-reflux valve and pour spout lid	50	

RECEPTAL® —20	000 Series (Liner capacity is approximate) - Latex-Free	
List No.	Description	Case Size
43445-01	2000 mL Receptal® eanister (for use with 2000 mL liners)	10
List No.	Description	Case Size
43027-01	2000 mL liner with pour spout lid	50
List No.	Description	Case Size
43041-01	2000 mL liner with non-pour spout (orthopedic liner) lid	50
List No.	Description	Case Size
43042-01	2000 mL liner with VAC-GARD® and non-pour spout lid	50
List No.	Description	Case Size
43044-05	2000 mL liner with VAC-GARD® and pour spout lid	50
List No.	Description	Case Size
43306-01	Receptal® Accu-Measure® device	12

EZE-VAC® Rigid Di	sposable Canister System	
List No.	Description	Case Size
43203-01	1200 mL canister with lid	32
List No.	Description	Case Size
43206-01	2000 mL canister with lid	24
List No.	Description	Case Size
43207-01	700 mL canister with lid	32
List No.	Description	Case Size
43213-01	3300 mL canister with lid (Tall, tapered and stackable)	12

RECEPTAL® & E	ZE-VAC $^{ ext{@}}$ Floorstands and Brackets. Most are stainless steel – Late	ex-Free
List No.	Description	Case Size
43444-05	Wall channel bracket for mounting canister (Receptal® or EZE-VAC®)	5
List No.	Description	Case Size
43474-01	Large head stainless steel floor stand (holds up to 12 Receptal® or 8 EZE-VAC® canisters)	1
List No.	Description	Case Size
43475-01	Small head stainless steel floor stand (holds up to 6 Receptal® or 4 EZE-VAC® canisters)	1
List No.	Description	Case Size
43506-02	Small head floor stand (holds up to 6 Receptal® or 4 EZE-VAC® canisters)	1
List No.	Description	Case Size
44303-10	Bracket clamp for adapting 43444-05 to stands, carts or IV poles	10
List No.	Description	Case Size
44304-05	Canister support bracket for mounting a regulator	5
List No.	Description	Case Size
44309-05	Regulator mount for canister support	5

RECEPTAL® & E	ZE-VAC® Accessories - Latex-Free	
List No.	Description	Case Size
43311-01	72" patient tubing	50
List No.	Description	Case Size
43315-25	3/8" pour spout adapter	25
List No.	Description	Case Size
43455-01	100' PVC suction tubing 0.3" I.D. (for canister to vacuum connection)	1
List No.	Description	Case Size
43457-25	24" PVC suction tubing 0.3" I.D. (for canister to vacuum connection)	25
Other Accessori	es	
List No.	Description	Case Size
43303-01	Specimen trap 4" - pour spout cover with 3/8" collection port, tethered cap and no stiffener	25
List No.	Description	Case Size
43303-02	Specimen trap 4" pour spout cover with 3/8" collection port, tethered cap and stiffener	25
List No.	Description	Case Size
43303-05	Specimen trap 8" - pour spout cover with 1/2" collection port, tethered cap and stiffener	25
List No.	Description	Case Size
43303-06	Specimen trap 8" - pour spout cover with 3/8" collection port, tethered cap and stiffener	25
List No.	Description	Case Size
43305-10	Canister attachment bracket for Accu-Measure® (Latex-Free)	10
List No.	Description	Case Size
43309-01	External VAC-GARD® filter/shutoff valve (Latex-Free)	100

List No. 43497-01	Description 3 way stopcock (for vacuum flow direction/control Latex-Free)	Case Size
List No. 43498-01	Description 2 way stopcock (for vacuum flow direction/control – Latex-Free)	Case Size
List No. 43501-05	Description 2/3 way stopcock holder (Latex-Free)	Case Size

negulator Moul	nts (Regulators not available from Hospira)	
List No. 43308-05	Description DISS Quick Connect Device (for use on 43423-01 for mounting a regulator with 43417-01)	Case Size 5
list No.	Description	Case Size
43400-05	Male DISS x 1/4" male Receptal® taper (for use with 43417–05 or when a regulator is wall mounted to vacuum and a simple taper is needed to attach tubing to regulator)	5
List No.	Description	Case Size
43404-05	Double 1/4" male Receptal® taper with 1 1/8" threaded shank and two locknuts (for use on floor stands to attach tubing from the floor stand to vacuum source)	5
List No.	Description	Case Size
43407-05	1/4" male pipe thread x $1/4$ " male Receptal® taper with 1 $1/8$ " threaded shank and two locknuts (for use for competitive regulators in the same way as $43408-05$)	5
List No.	Description	Case Size
43408-05	1/8" male pipe thread x 1/4" male Receptal® taper with 1 1/8" threaded shank and two locknuts	5

Regulator Mounts (Regulators not available from Hospira) (cont.)				
List No.	Description	Case Size		
43417-05	Female DISS \times 1/8" male pipe thread (for use to insert into the bottom of a regulator and is then attached to 43401-05 for regulator mounting on a floor stand or wall bracket)	5		
List No.	Description	Case Size		
43420-05	This is the same as list #43417-01, except it is a 1/4" male pipe thread and would be used for mounting competitive regulators that are 1/4" in size with list #43401-01	5		
List No.	Description	Case Size		
43421-05	1/4" male pipe thread x 1/4" male Receptal® taper	5		
List No.	Description	Case Size		
43422-05	1/8" male pipe thread x 1/4" male Receptal® taper	5		
List No.	Description	Case Size		
43491-05	Male DISS \times 1/8" NPT adapter (for fitting a regulator guage with DISS capability. For use with 43203-01 and 43201-01)	5		