

INNOVATIVE PRODUCTS FOR G, J AND GJ TUBES



AVANOS Medical, Your Preferred Enteral Partner, offers a wide variety of innovative, high-quality enteral feeding tubes and accessories uniquely designed for delivering nutrition and medication to pediatric and adult patients.

As the pioneer in the development of the first balloon retained tube designed for gastrostomy feeding, AVANOS continues to manufacture tubes and accessories with the highest standards of workmanship at our state-of-the-art manufacturing facilities. Our products provide consistent performance. Our products are designed for safety and convenience. We stand behind our products with knowledgeable sales and customer service professionals. Today, our MIC*/MIC-KEY* brand of enteral feeding tubes is well recognized and accepted among healthcare professionals for intermediate and long-term feeding applications.

MIC*/MIC-KEY*

Enteral Feeding Tubes & Accessories

INNOVATIVE PRODUCT DESIGNS


MIC*/MIC-KEY* Enteral Feeding Tubes and Accessories include a wide variety of innovative product designs for gastric, jejunal and gastro-jejunal access: from the lock-and-key design of the MIC-KEY* Feeding Tube extension set connection mechanism, to the one-piece MIC-KEY* GJ Feeding Tube that allows simultaneous gastric decompression and jejunal feeding.

MOST MIC*/MIC-KEY* ENTERAL FEEDING TUBES INCLUDE:

- High-clarity, medical-grade silicone construction, designed for visibility and drapability
- Ventilated SECUR-LOK* ring, the external retention mechanism design allows air to circulate around the stoma site and reduce tube pressure 360° from the sides of the stoma wall
- A patented universal connector, designed to preserve tube life and minimize unintentional disconnects
- An inflatable internal retention balloon, which reduces patient trauma experienced with obturated devices

THE RIGHT CHOICE

Helping your patients face the use of a gastrostomy tube, possibly for a lifetime, can be a challenge. AVANOS's variety of quality feeding tubes and accessories helps you make the right choice for your patient. By selecting MIC*/MIC-KEY* Enteral Feeding Tubes, you may be assured that you have chosen innovative products that consistently perform and are well recognized and accepted among healthcare professionals.



Patented Universal Connector

SECUR-LOK* External Retention Ring

Inflatable Internal Retention Balloon

GASTRIC ACCESS FEEDING TUBES

AVANOS[®] Gastric Access Feeding Tubes are designed for the delivery of enteral nutrition to the stomach and are available in low-profile and conventional, standard-length designs.




The **MIC-KEY[®] G FEEDING TUBE** is a step up from conventional gastrostomy feeding tubes. Highly popular and widely prescribed, the MIC-KEY[®] G Feeding Tube is unobtrusive and easy to conceal, making it an ideal feeding tube for virtually all individuals. Specific features include:

- Medical-Grade Silicone Construction
- Low-Profile Design
- Tapered Distal Tip
- Inflatable Silicone Internal Retention Balloon
- Distal Tip Recessed at 5ml
- Proximal Anti-Reflux Valve
- SECUR-LOK[®] Extension Set Connector Mechanism
- Radiopaque Stripe



The first of its kind, the **MIC[®] G FEEDING TUBE** is indicated for patients requiring continuous enteral feeding. Specific features include:

- Medical-Grade Silicone Construction
- Recessed Distal Tip at Recommended Fill Volume
- Inflatable Silicone Internal Retention Balloon
- SECUR-LOK[®] External Retention Ring
- Universal Feeding Port Connector
- Medication Port
- Tapered Distal Tip
- Dual Exit Ports
- Radiopaque Stripe



The **MIC[®] BOLUS G FEEDING TUBE** is utilized for patients who are on a bolus feeding regimen. Specific features include:

- Medical-Grade Silicone Construction
- Recessed Distal Tip at Recommended Fill Volume
- Inflatable Silicone Internal Retention Balloon
- SECUR-LOK[®] External Retention Ring
- Tapered Distal Tip
- Dual Exit Ports
- Radiopaque Stripe




The **MIC[®] PEG FEEDING TUBE KIT** is designed for placement convenience and comes in complete placement kits for both the PULL and PUSH placement methods. Available in 14, 20 and 24 French sizes. Specific features include:

- Medical-Grade Silicone Construction
- QUAD-LOK[®] External Retention Ring
- Universal and Bolus Feeding Port Connectors
- Medication Port
- Collapsible Internal Retention Bumper
- Radiopaque Stripe and Bumper
- Tubing Clamp
- Available in PUSH (OTW) and PULL Placement Methods


JEJUNAL ACCESS FEEDING TUBES

AVANOS® Jejunol Access Feeding Tubes are designed for the delivery of enteral nutrition to the small bowel and are available in low-profile and conventional designs.




The **MIC-KEY® J FEEDING TUBE** is indicated for patients requiring jejunal feeding only, but who want an inconspicuous external segment. With an external bolster designed like the widely prescribed MIC-KEY® G Feeding Tube, specific features include:

- Medical-Grade Silicone Construction
- Low-Profile Design
- Tapered Distal Tip
- Inflatable Silicone Internal Retention Balloon
- Distal Tip Recessed at 5ml
- Proximal Anti-Reflux Valve
- SECUR-LOK® Extension Set Connector Mechanism
- Trimmable Distal Tip
- Radiopaque Stripe
- Wide Variety of Extension Sets Available



The **MIC® J FEEDING TUBE** is indicated for patients requiring jejunal feeding only. The distal tip can be trimmed to suit individual patient needs. Specific features include:

- Medical-Grade Silicone Construction
- Inflatable Silicone Internal Retention Balloon
- SECUR-LOK® External Retention Ring
- Universal Feeding Port Connector
- Medication Port
- Trimmable Distal Tip
- Radiopaque Stripe
- Gamma Sterilized



The **MIC® JEJUNOSTOMY FEEDING TUBE** is indicated for patients requiring jejunal feeding. It is surgically placed directly into the small bowel, utilizing the Witzel Tunnel technique for secure placement and to help minimize leakage. Featuring a large 9 French internal diameter to help prevent clogging, specific features include:

- Medical-Grade Silicone Construction
- Trimmable Distal Tip
- Suture Wings
- Large 9 French Internal Diameter
- Radiopaque Stripe


GASTRIC-JEJUNAL ACCESS FEEDING TUBES

AVANOS[®] Gastro-Jejunal Access Feeding Tubes are designed for the delivery of enteral nutrition to the small bowel while allowing simultaneous decompression of the stomach.




The unique **MIC-KEY[®] GJ FEEDING TUBE** is also designed for patients who require simultaneous gastric decompression and jejunal feeding. Available in several sizes to meet individual patient needs. Specific features include:

- Medical-Grade Silicone Construction
- Inflatable Silicone Internal Retention Balloon
- Universal Jejunal Feeding Port Connector
- Multiple Gastric & Jejunal Exit Ports
- Gastric Decompression Port
- Tungsten Weighted Jejunal Portion
- Tapered Distal Tip
- Radiopaque Stripe



The unique **MIC[®] GASTRIC-JEJUNAL FEEDING TUBE** is also designed for patients who require simultaneous gastric decompression and jejunal feeding. Available in several sizes to meet individual patient needs. Specific features include:

- Medical-Grade Silicone Construction
- Inflatable Silicone Internal Retention Balloon
- SECUR-LOK[®] External Retention Ring
- Universal Jejunal Feeding Port Connector
- Gastric Decompression Port
- Multiple Gastric & Jejunal Exit Ports
- Tungsten Weighted Jejunal Portion
- Tapered Distal Tip
- Radiopaque Stripe



The **MIC[®] GASTRO-ENTERIC FEEDING TUBE** is a precursor to the very popular MIC[®] GJ Feeding Tube. Also designed for simultaneous gastric decompression and jejunal feeding, it is indicated for patients where smaller distal jejunal segment may be indicated. Specific features include:

- Medical-Grade Silicone Construction
- Inflatable Silicone Internal Retention Balloon
- SECUR-LOK[®] External Retention Ring
- Universal Jejunal Feeding Port Connector
- Gastric Decompression Port
- Multiple Gastric & Jejunal Exit Ports
- Radiopaque Stripe

ACCESSORIES

MIC[®]/MIC-KEY[®] Accessories support gastric, jejunal and gastro-jejunal access feeding tubes.



- SECUR-LOK[®] "Key" Connection Mechanism
- Medication Port on Non-Bolus Extension Sets
- Individually Packaged Extension Sets
- Tubing Clamp

MIC*/MIC-KEY*

Enteral Feeding Tubes & Accessories

ORDERING CODES

GASTRIC ACCESS FEEDING TUBES

AVANOS* MIC-KEY* G FEEDING TUBES/KITS

Legacy Stock #	ENFit® Stock #	Size (Outer Dia.)	Stoma Length	Balloon Volume	Packaging
0120-12-0.8	8140-12-0.8	12 Fr	0.8cm	3ml Balloon	1 Each
0120-12-1.0	8140-12-1.0	12 Fr	1.0cm	3ml Balloon	1 Each
0120-12-1.2	8140-12-1.2	12 Fr	1.2cm	3ml Balloon	1 Each
0120-12-1.5	8140-12-1.5	12 Fr	1.5cm	3ml Balloon	1 Each
0120-12-1.7	8140-12-1.7	12 Fr	1.7cm	3ml Balloon	1 Each
0120-12-2.0	8140-12-2.0	12 Fr	2.0cm	3ml Balloon	1 Each
0120-12-2.3	8140-12-2.3	12 Fr	2.3cm	3ml Balloon	1 Each
0120-12-2.5	8140-12-2.5	12 Fr	2.5cm	3ml Balloon	1 Each
0120-12-2.7	8140-12-2.7	12 Fr	2.7cm	3ml Balloon	1 Each
0120-12-3.0	8140-12-3.0	12 Fr	3.0cm	3ml Balloon	1 Each
0120-12-3.5	8140-12-3.5	12 Fr	3.5cm	3ml Balloon	1 Each
0120-12-4.0	8140-12-4.0	12 Fr	4.0cm	3ml Balloon	1 Each
0120-14-0.8	8140-14-0.8	14 Fr	0.8cm	5ml Balloon	1 Each
0120-14-1.0	8140-14-1.0	14 Fr	1.0cm	5ml Balloon	1 Each
0120-14-1.2	8140-14-1.2	14 Fr	1.2cm	5ml Balloon	1 Each
0120-14-1.5	8140-14-1.5	14 Fr	1.5cm	5ml Balloon	1 Each
0120-14-1.7	8140-14-1.7	14 Fr	1.7cm	5ml Balloon	1 Each
0120-14-2.0	8140-14-2.0	14 Fr	2.0cm	5ml Balloon	1 Each
0120-14-2.3	8140-14-2.3	14 Fr	2.3cm	5ml Balloon	1 Each
0120-14-2.5	8140-14-2.5	14 Fr	2.5cm	5ml Balloon	1 Each
0120-14-2.7	8140-14-2.7	14 Fr	2.7cm	5ml Balloon	1 Each
0120-14-3.0	8140-14-3.0	14 Fr	3.0cm	5ml Balloon	1 Each
0120-14-3.5	8140-14-3.5	14 Fr	3.5cm	5ml Balloon	1 Each
0120-14-4.0	8140-14-4.0	14 Fr	4.0cm	5ml Balloon	1 Each
0120-14-4.5	8140-14-4.5	14 Fr	4.5cm	5ml Balloon	1 Each
0120-14-5.0	8140-14-5.0	14 Fr	5.0cm	5ml Balloon	1 Each
0120-16-0.8	8140-16-0.8	16 Fr	0.8cm	5ml Balloon	1 Each
0120-16-1.0	8140-16-1.0	16 Fr	1.0cm	5ml Balloon	1 Each
0120-16-1.2	8140-16-1.2	16 Fr	1.2cm	5ml Balloon	1 Each
0120-16-1.5	8140-16-1.5	16 Fr	1.5cm	5ml Balloon	1 Each
0120-16-1.7	8140-16-1.7	16 Fr	1.7cm	5ml Balloon	1 Each
0120-16-2.0	8140-16-2.0	16 Fr	2.0cm	5ml Balloon	1 Each
0120-16-2.3	8140-16-2.3	16 Fr	2.3cm	5ml Balloon	1 Each
0120-16-2.5	8140-16-2.5	16 Fr	2.5cm	5ml Balloon	1 Each
0120-16-2.7	8140-16-2.7	16 Fr	2.7cm	5ml Balloon	1 Each
0120-16-3.0	8140-16-3.0	16 Fr	3.0cm	5ml Balloon	1 Each
0120-16-3.5	8140-16-3.5	16 Fr	3.5cm	5ml Balloon	1 Each
0120-16-4.0	8140-16-4.0	16 Fr	4.0cm	5ml Balloon	1 Each

Legacy Stock #	ENFit® Stock #	Size (Outer Dia.)	Stoma Length	Balloon Volume	Packaging
0120-16-4.5	8140-16-4.5	16 Fr	4.5cm	5ml Balloon	1 Each
0120-16-5.0	8140-16-5.0	16 Fr	5.0cm	5ml Balloon	1 Each
0120-18-0.8	8140-18-0.8	18 Fr	0.8cm	5ml Balloon	1 Each
0120-18-1.0	8140-18-1.0	18 Fr	1.0cm	5ml Balloon	1 Each
0120-18-1.2	8140-18-1.2	18 Fr	1.2cm	5ml Balloon	1 Each
0120-18-1.5	8140-18-1.5	18 Fr	1.5cm	5ml Balloon	1 Each
0120-18-1.7	8140-18-1.7	18 Fr	1.7cm	5ml Balloon	1 Each
0120-18-2.0	8140-18-2.0	18 Fr	2.0cm	5ml Balloon	1 Each
0120-18-2.3	8140-18-2.3	18 Fr	2.3cm	5ml Balloon	1 Each
0120-18-2.5	8140-18-2.5	18 Fr	2.5cm	5ml Balloon	1 Each
0120-18-2.7	8140-18-2.7	18 Fr	2.7cm	5ml Balloon	1 Each
0120-18-3.0	8140-18-3.0	18 Fr	3.0cm	5ml Balloon	1 Each
0120-18-3.5	8140-18-3.5	18 Fr	3.5cm	5ml Balloon	1 Each
0120-18-4.0	8140-18-4.0	18 Fr	4.0cm	5ml Balloon	1 Each
0120-18-4.5	8140-18-4.5	18 Fr	4.5cm	5ml Balloon	1 Each
0120-18-5.0	8140-18-5.0	18 Fr	5.0cm	5ml Balloon	1 Each
0120-20-0.8	8140-20-0.8	20 Fr	0.8cm	5ml Balloon	1 Each
0120-20-1.0	8140-20-1.0	20 Fr	1.0cm	5ml Balloon	1 Each
0120-20-1.2	8140-20-1.2	20 Fr	1.2cm	5ml Balloon	1 Each
0120-20-1.5	8140-20-1.5	20 Fr	1.5cm	5ml Balloon	1 Each
0120-20-1.7	8140-20-1.7	20 Fr	1.7cm	5ml Balloon	1 Each
0120-20-2.0	8140-20-2.0	20 Fr	2.0cm	5ml Balloon	1 Each
0120-20-2.3	8140-20-2.3	20 Fr	2.3cm	5ml Balloon	1 Each
0120-20-2.5	8140-20-2.5	20 Fr	2.5cm	5ml Balloon	1 Each
0120-20-2.7	8140-20-2.7	20 Fr	2.7cm	5ml Balloon	1 Each
0120-20-3.0	8140-20-3.0	20 Fr	3.0cm	5ml Balloon	1 Each
0120-20-3.5	8140-20-3.5	20 Fr	3.5cm	5ml Balloon	1 Each
0120-20-4.0	8140-20-4.0	20 Fr	4.0cm	5ml Balloon	1 Each
0120-20-4.5	8140-20-4.5	20 Fr	4.5cm	5ml Balloon	1 Each
0120-20-5.0	8140-20-5.0	20 Fr	5.0cm	5ml Balloon	1 Each
0120-24-1.5	8140-24-1.5	24 Fr	1.5cm	5ml Balloon	1 Each
0120-24-1.7	8140-24-1.7	24 Fr	1.7cm	5ml Balloon	1 Each
0120-24-2.0	8140-24-2.0	24 Fr	2.0cm	5ml Balloon	1 Each
0120-24-2.3	8140-24-2.3	24 Fr	2.3cm	5ml Balloon	1 Each
0120-24-2.5	8140-24-2.5	24 Fr	2.5cm	5ml Balloon	1 Each
0120-24-2.7	8140-24-2.7	24 Fr	2.7cm	5ml Balloon	1 Each
0120-24-3.0	8140-24-3.0	24 Fr	3.0cm	5ml Balloon	1 Each
0120-24-3.5	8140-24-3.5	24 Fr	3.5cm	5ml Balloon	1 Each
0120-24-4.0	8140-24-4.0	24 Fr	4.0cm	5ml Balloon	1 Each
0120-24-4.5	8140-24-4.5	24 Fr	4.5cm	5ml Balloon	1 Each
0120-24-5.0	8140-24-5.0	24 Fr	5.0cm	5ml Balloon	1 Each

AVANOS* MIC* G FEEDING TUBES

Legacy Stock #	Size (Outer Dia.)	Balloon Volume	Packaging
0100-12LV†	12 Fr	3-5ml Balloon	1 Each
0100-14LV	14 Fr	3-5ml Balloon	1 Each
0100-16LV	16 Fr	3-5ml Balloon	1 Each
0100-18	18 Fr	7-10ml Balloon	1 Each
0100-20	20 Fr	7-10ml Balloon	1 Each
0100-22	22 Fr	7-10ml Balloon	1 Each
0100-24	24 Fr	7-10ml Balloon	1 Each
0100-26	26 Fr	7-10ml Balloon	1 Each
0100-28	28 Fr	7-10ml Balloon	1 Each
0100-30	30 Fr	7-10ml Balloon	1 Each

AVANOS* MIC* BOLUS G FEEDING TUBES

Legacy Stock #	Size (Outer Dia.)	Balloon Volume	Packaging
0110-12LV†	12 Fr	3-5ml Balloon	1 Each
0110-14LV	14 Fr	3-5ml Balloon	1 Each
0110-16LV	16 Fr	3-5ml Balloon	1 Each
0110-18	18 Fr	7-10ml Balloon	1 Each
0110-20	20 Fr	7-10ml Balloon	1 Each
0110-22	22 Fr	7-10ml Balloon	1 Each
0110-24	24 Fr	7-10ml Balloon	1 Each

AVANOS* MIC* PERCUTANEOUS ENDOSCOPIC GASTROSTOMY (PEG) FEEDING TUBES/KITS

PUSH PEG Kit Includes:

- 1- Traction Removable PUSH PEG Tube
- 1- 5 Fr Retrieval Snare (1.8mm X 260cm)
- 1- Universal Feeding Adapter
- 1- External Retention Bolster
- 1- 5ml Vial 1% Lidocaine HCl
- 1- 19ga X 1.5in. Filter Needle
- 1- 25ga X 1.5in. Needle
- 1- Povidone Iodine Ointment
- 1- Fenestrated Drape with Adhesive Tabs
- 1- 6ml Syringe
- 1- Scalpel with #11 Blade
- 1- Stainless Steel Introducer Needle
- 1- Guidewire
- 1- Stainless Steel Scissors
- 1- Stainless Steel Hemostat
- 1- Tubing Clamp
- 1- Water-Soluble Lubricant
- 4- 4in. X 4in. Gauze Sponges
- 2- 2in. X 2in. Slit Gauze Sponges
- 1- Chloraprep® Swabstick Packet
- 1- Lidocaine Label

PULL PEG Kit Includes:

- 1- Traction Removable PULL PEG Tube
- 1- 5 Fr Retrieval Snare (1.8mm X 260cm)
- 1- Universal Feeding Adapter
- 1- External Retention Bolster
- 1- 5ml Vial 1% Lidocaine HCl
- 1- 19ga X 1.5in. Filter Needle
- 1- 25ga X 1.5in. Needle
- 1- Povidone Iodine Swabstick Packet
- 1- Povidone Iodine Ointment
- 1- Fenestrated Drape with Adhesive Tabs
- 1- 6ml Syringe
- 1- Scalpel with #11 Blade
- 1- Introducer Needle
- 1- Looped Placement Wire
- 1- Stainless Steel Scissors
- 1- Stainless Steel Hemostat
- 1- Tubing Clamp
- 1- Water-Soluble Lubricant
- 4- 4in. X 4in. Gauze Sponges
- 2- 2in. X 2in. Slit Gauze Sponges
- 1- Chloraprep® Swabstick Packet
- 1- Lidocaine Label

Legacy Stock #	ENFit® Stock #	Size (Outer Dia.)	Balloon Volume	Packaging
7150-14	8150-14	14 Fr	PUSH (Over-the-Wire) Method	2 Each
7150-20	8150-20	20 Fr	PUSH (Over-the-Wire) Method	2 Each
7150-24	8150-24	24 Fr	PUSH (Over-the-Wire) Method	2 Each
7160-14	8160-14	14 Fr	PULL Method	2 Each
7160-20	8160-20	20 Fr	PULL Method	2 Each
7160-24	8160-24	24 Fr	PULL Method	2 Each

†Non-Recessed Tip

AVANOS* MIC* SAFETY PEG FEEDING TUBES/KITS

PUSH PEG Kit Includes:

- 1- Traction Removable PUSH PEG Tube
- 1- 5 Fr Retrieval Snare (1.8mm X 260cm)
- 1- Universal Feeding Adapter
- 1- External Retention Bolster
- 1- 5ml Vial 1% Lidocaine HCl
- 1- 18ga X 1.5in. 5 Micron Blunt Filter Needle
- 1- 22ga X 1.5in. Needle
- 1- Chloraprep® Swabstick Packet
- 1- Povidone Iodine Ointment
- 1- Fenestrated Drape with Adhesive Tabs
- 1- Syringe, 6ml
- 1- Safety Scalpel with #11 Blade
- 1- Introducer Safety Needle, 17 gauge
- 1- Guidewire (.038in. X 260cm)
- 1- Stainless Steel Scissors
- 1- Stainless Steel Hemostat
- 1- Tubing Clamp
- 1- Water-Soluble Lubricant
- 4- 4in. X 4in. Gauze Sponges
- 2- 2in. X 2in. Slit Gauze Sponges
- 1- Lidocaine Label

PULL PEG Kit Includes:

- 1- Traction Removable PULL PEG Tube
- 1- 5 Fr Retrieval Snare (1.8mm X 260cm)
- 1- Universal Feeding Adapter
- 1- External Retention Bolster
- 1- 5ml Vial 1% Lidocaine HCl
- 1- 18ga X 1.5in. 5 Micron Blunt Filter Needle
- 1- 22ga X 1.5in. Needle
- 1- Chloraprep® Swabstick Packet
- 1- Povidone Iodine Ointment
- 1- Fenestrated Drape with Adhesive Tabs
- 1- Syringe, 6ml
- 1- Safety Scalpel with #11 Blade
- 1- Introducer Safety Needle, 17 gauge
- 1- Looped Placement Wire
- 1- Stainless Steel Scissors
- 1- Stainless Steel Hemostat
- 1- Tubing Clamp
- 1- Water-Soluble Lubricant
- 4- 4in. X 4in. Gauze Sponges
- 2- 2in. X 2in. Slit Gauze Sponges
- 1- Lidocaine Label

Legacy Stock #	ENFit® Stock #	Size (Outer Dia)	Balloon Volume	Packaging
7170-14	8170-14	14 Fr	PUSH (Over-the-Wire) Method	2 Each
7170-20	8170-20	20 Fr	PUSH (Over-the-Wire) Method	2 Each
7170-24	8170-24	24 Fr	PUSH (Over-the-Wire) Method	2 Each
7180-14	8180-14	14 Fr	PULL Method	2 Each
7180-20	8180-20	20 Fr	PULL Method	2 Each
7180-24	8180-24	24 Fr	PULL Method	2 Each

JEJUNAL ACCESS FEEDING TUBES

AVANOS* MIC-KEY* J FEEDING TUBES/KITS

Legacy Stock #	ENFit® Stock #	Size (Outer Dia.)	Stoma Length	Balloon Volume	Packaging
0230-14-0.8	8230-14-0.8	14 Fr	0.8cm	5ml Balloon	1 Each
0230-14-1.0	8230-14-1.0	14 Fr	1.0cm	5ml Balloon	1 Each
0230-14-1.2	8230-14-1.2	14 Fr	1.2cm	5ml Balloon	1 Each
0230-14-1.5	8230-14-1.5	14 Fr	1.5cm	5ml Balloon	1 Each
0230-14-1.7	8230-14-1.7	14 Fr	1.7cm	5ml Balloon	1 Each
0230-14-2.0	8230-14-2.0	14 Fr	2.0cm	5ml Balloon	1 Each
0230-14-2.3	8230-14-2.3	14 Fr	2.3cm	5ml Balloon	1 Each
0230-14-2.5	8230-14-2.5	14 Fr	2.5cm	5ml Balloon	1 Each
0230-14-2.7	8230-14-2.7	14 Fr	2.7cm	5ml Balloon	1 Each
0230-14-3.0	8230-14-3.0	14 Fr	3.0cm	5ml Balloon	1 Each
0230-14-3.5	8230-14-3.5	14 Fr	3.5cm	5ml Balloon	1 Each
0230-14-4.0	8230-14-4.0	14 Fr	4.0cm	5ml Balloon	1 Each
0230-14-4.5	8230-14-4.5	14 Fr	4.5cm	5ml Balloon	1 Each
0230-18-0.8	8230-18-0.8	18 Fr	0.8cm	5ml Balloon	1 Each
0230-18-1.0	8230-18-1.0	18 Fr	1.0cm	5ml Balloon	1 Each
0230-18-1.2	8230-18-1.2	18 Fr	1.2cm	5ml Balloon	1 Each
0230-18-1.5	8230-18-1.5	18 Fr	1.5cm	5ml Balloon	1 Each
0230-18-1.7	8230-18-1.7	18 Fr	1.7cm	5ml Balloon	1 Each
0230-18-2.0	8230-18-2.0	18 Fr	2.0cm	5ml Balloon	1 Each
0230-18-2.3	8230-18-2.3	18 Fr	2.3cm	5ml Balloon	1 Each
0230-18-2.5	8230-18-2.5	18 Fr	2.5cm	5ml Balloon	1 Each
0230-18-2.7	8230-18-2.7	18 Fr	2.7cm	5ml Balloon	1 Each
0230-18-3.0	8230-18-3.0	18 Fr	3.0cm	5ml Balloon	1 Each
0230-18-3.5	8230-18-3.5	18 Fr	3.5cm	5ml Balloon	1 Each
0230-18-4.0	8230-18-4.0	18 Fr	4.0cm	5ml Balloon	1 Each
0230-18-4.5	8230-18-4.5	18 Fr	4.5cm	5ml Balloon	1 Each

AVANOS* MIC* J FEEDING TUBES

Legacy Stock #	Size (Outer Dia.)	Balloon Volume	Packaging
0200-12LV	12 Fr	3-5ml Balloon	1 Each
0200-14	14 Fr	7-10ml Balloon	1 Each
0200-16	16 Fr	7-10ml Balloon	1 Each
0200-18	18 Fr	7-10ml Balloon	1 Each
0200-20	20 Fr	7-10ml Balloon	1 Each
0200-22	22 Fr	7-10ml Balloon	1 Each
0200-24	24 Fr	7-10ml Balloon	1 Each

AVANOS* MIC* JEJUNOSTOMY FEEDING TUBE

Includes:

- 1- Jejunostomy Tube
- 1- Subcutaneous Extraducer
- 1- Feeding Port Adapter
- 1- Luer-Lok® Cap

Legacy Stock #	Outer Dia./ Inner Dia.	Packaging
0301-14	14 Fr/9 Fr	1 Each

GASTRIC-JEJUNAL ACCESS FEEDING TUBES

AVANOS* MIC* GJ FEEDING TUBES – ENDOSCOPIC/RADIOLOGY PLACEMENT KITS

Kit Includes:

- 1- GJ Feeding Tube
- 1- Silicone-Coated Guidewire
- 1- Luer-Slip™ Syringe – Prefilled with 10ml water
- 1- 35ml Catheter Tip Syringe
- 1- Water-Soluble Lubricant
- 1- SECUR-LOK® External Retention Ring
- 1- Directions for Use
- 1- Patient Care Guidelines

Legacy Stock #	Size (Outer Dia.)	Jejunal Length+	Balloon Volume	Packaging
0250-16-15	16 Fr	15cm	7-10ml Balloon	1 Each
0250-16-22	16 Fr	22cm	7-10ml Balloon	1 Each
0250-16-30	16 Fr	30cm	7-10ml Balloon	1 Each
0250-16	16 Fr	45cm	7-10ml Balloon	1 Each
0250-18-22	18 Fr	22cm	7-10ml Balloon	1 Each
0250-18-30	18 Fr	30cm	7-10ml Balloon	1 Each
0250-18	18 Fr	45cm	7-10ml Balloon	1 Each
0250-22	22 Fr	45cm	7-10ml Balloon	1 Each

AVANOS* MIC* GJ FEEDING TUBES – SURGICAL PLACEMENT KITS

Kit Includes:

- 1- GJ Feeding Tube
- 1- Abdominal Trocar and Sheath
- 1- Slotted Transpyloric Cannula
- 1- Luer-Slip™ Syringe – Prefilled with 10ml water
- 1- 35ml Catheter Tip Syringe
- 1- Water-Soluble Lubricant
- 1- 2in. X 2in. Slit Gauze Sponge
- 1- SECUR-LOK® External Retention Ring
- 1- Directions for Use
- 1- Postoperative Care Guidelines
- 1- Patient Care Guidelines

Legacy Stock #	Size (Outer Dia.)	Jejunal Length+	Balloon Volume	Packaging
0260-16	16 Fr	45cm	7-10ml Balloon	1 Each
0260-18	18 Fr	45cm	7-10ml Balloon	1 Each
0260-22	22 Fr	45cm	7-10ml Balloon	1 Each

AVANOS* MIC-KEY* GJ FEEDING TUBES – ENDOSCOPIC/RADIOLOGY PLACEMENT KITS

Kit Includes:

- 1- MIC-KEY® Low-Profile GJ Feeding Tube
- 1- Introducer Cannula
- 1- MIC-KEY® Extension Set with SECUR-LOK® Right-Angle Connector and 2 Port “Y” and Clamp – 12in. Length
- 1- MIC-KEY® Bolus Extension Set with Catheter Tip, SECUR-LOK® Straight Connector and Clamp – 12in. Length
- 1- 6ml Luer-Slip™ Syringe
- 1- 35ml Catheter Tip Syringe
- 4- Gauze Pads
- 1- Directions For Use
- 1- Patient Care Guidelines

AVANOS* MIC-KEY* GJ FEEDING TUBES – ENDOSCOPIC/RADIOLOGY PLACEMENT KITS

Legacy Stock #	ENFit® Stock #	Size (Outer Dia.)	Stoma Length	Jejunal Length+	Balloon Volume	Packaging
0270-16-1.0-15	8270-16-1.0-15	16 Fr	1.0cm	15cm	3-5ml	1 Each
0270-16-1.0-22	8270-16-1.0-22	16 Fr	1.0cm	22cm	3-5ml	1 Each
0270-16-1.2-15	8270-16-1.2-15	16 Fr	1.2cm	15cm	3-5ml	1 Each
0270-16-1.2-22	8270-16-1.2-22	16 Fr	1.2cm	22cm	3-5ml	1 Each
0270-16-1.2-30	8270-16-1.2-30	16 Fr	1.2cm	30cm	3-5ml	1 Each
0270-16-1.5-15	8270-16-1.5-15	16 Fr	1.5cm	15cm	3-5ml	1 Each
0270-16-1.5-22	8270-16-1.5-22	16 Fr	1.5cm	22cm	3-5ml	1 Each
0270-16-1.5-30	8270-16-1.5-30	16 Fr	1.5cm	30cm	3-5ml	1 Each
0270-16-1.5-45	8270-16-1.5-45	16 Fr	1.5cm	45cm	7-10ml	1 Each
0270-16-1.7-15	8270-16-1.7-15	16 Fr	1.7cm	15cm	3-5ml	1 Each
0270-16-1.7-22	8270-16-1.7-22	16 Fr	1.7cm	22cm	3-5ml	1 Each
0270-16-1.7-30	8270-16-1.7-30	16 Fr	1.7cm	30cm	3-5ml	1 Each
0270-16-1.7-45	8270-16-1.7-45	16 Fr	1.7cm	45cm	7-10ml	1 Each
0270-16-2.0-22	8270-16-2.0-22	16 Fr	2.0cm	22cm	3-5ml	1 Each
0270-16-2.0-30	8270-16-2.0-30	16 Fr	2.0cm	30cm	3-5ml	1 Each
0270-16-2.0-45	8270-16-2.0-45	16 Fr	2.0cm	45cm	7-10ml	1 Each
0270-16-2.3-30	8270-16-2.3-30	16 Fr	2.3cm	30cm	3-5ml	1 Each
0270-16-2.3-45	8270-16-2.3-45	16 Fr	2.3cm	45cm	7-10ml	1 Each
0270-16-2.5-30	8270-16-2.5-30	16 Fr	2.5cm	30cm	3-5ml	1 Each
0270-16-2.5-45	8270-16-2.5-45	16 Fr	2.5cm	45cm	7-10ml	1 Each
0270-16-2.7-30	8270-16-2.5-45	16 Fr	2.7cm	30cm	3-5ml	1 Each
0270-16-2.7-45	8270-16-2.7-45	16 Fr	2.7cm	45cm	7-10ml	1 Each
0270-16-3.0-45	8270-16-3.0-45	16 Fr	3.0cm	45cm	7-10ml	1 Each
0270-16-3.5-45	8270-16-3.5-45	16 Fr	3.5cm	45cm	7-10ml	1 Each
0270-18-1.2-22	8270-18-1.2-22	18 Fr	1.2cm	22cm	3-5ml	1 Each
0270-18-1.2-30	8270-18-1.2-30	18 Fr	1.2cm	30cm	3-5ml	1 Each
0270-18-1.2-22	8270-18-1.2-22	18 Fr	1.2cm	22cm	3-5ml	1 Each
0270-18-1.2-30	8270-18-1.2-30	18 Fr	1.2cm	30cm	3-5ml	1 Each
0270-18-1.5-22	8270-18-1.5-22	18 Fr	1.5cm	22cm	3-5ml	1 Each
0270-18-1.5-30	8270-18-1.5-30	18 Fr	1.5cm	30cm	3-5ml	1 Each
0270-18-1.5-45	8270-18-1.5-45	18 Fr	1.5cm	45cm	7-10ml	1 Each
0270-18-1.7-22	8270-18-1.7-22	18 Fr	1.7cm	22cm	3-5ml	1 Each
0270-18-1.7-30	8270-18-1.7-30	18 Fr	1.7cm	30cm	3-5ml	1 Each
0270-18-1.7-45	8270-18-1.7-45	18 Fr	1.7cm	45cm	7-10ml	1 Each
0270-18-2.0-22	8270-18-2.0-22	18 Fr	2.0cm	22cm	3-5ml	1 Each
0270-18-2.0-30	8270-18-2.0-30	18 Fr	2.0cm	30cm	3-5ml	1 Each
0270-18-2.0-45	8270-18-2.0-45	18 Fr	2.0cm	45cm	7-10ml	1 Each
0270-18-2.3-30	8270-18-2.3-30	18 Fr	2.3cm	30cm	3-5ml	1 Each
0270-18-2.3-45	8270-18-2.3-45	18 Fr	2.3cm	45cm	7-10ml	1 Each
0270-18-2.5-30	8270-18-2.5-30	18 Fr	2.5cm	30cm	3-5ml	1 Each
0270-18-2.5-45	8270-18-2.5-45	18 Fr	2.5cm	45cm	7-10ml	1 Each
0270-18-2.7-30	8270-18-2.7-30	18 Fr	2.7cm	30cm	3-5ml	1 Each
0270-18-2.7-45	8270-18-2.7-45	18 Fr	2.7cm	45cm	7-10ml	1 Each
0270-18-3.0-45	8270-18-3.0-45	18 Fr	3.0cm	45cm	7-10ml	1 Each
0270-18-3.5-45	8270-18-3.5-45	18 Fr	3.5cm	45cm	7-10ml	1 Each

+Jejunal length measured from the distal aspect of the retention balloon to the catheter tip.

AVANOS* MIC-KEY* GJ FEEDING TUBES – ENDOSCOPIC/RADIOLOGY PLACEMENT KITS (CONT'D)

Legacy Stock #	ENFit® Stock #	Size (Outer Dia.)	Stoma Length	Jejunal Length+	Balloon Volume	Packaging
0270-22-1.5-45	8270-22-1.5-45	22 Fr	1.5cm	45cm	7-10ml	1 Each
0270-22-1.7-45	8270-22-1.7-45	22 Fr	1.7cm	45cm	7-10ml	1 Each
0270-22-2.0-45	8270-22-2.0-45	22 Fr	2.0cm	45cm	7-10ml	1 Each
0270-22-2.3-45	8270-22-2.3-45	22 Fr	2.3cm	45cm	7-10ml	1 Each
0270-22-2.5-45	8270-22-2.5-45	22 Fr	2.5cm	45cm	7-10ml	1 Each
0270-22-2.7-45	8270-22-2.7-45	22 Fr	2.7cm	45cm	7-10ml	1 Each
0270-22-3.0-45	8270-22-3.0-45	22 Fr	3.0cm	45cm	7-10ml	1 Each
0270-22-3.5-45	8270-22-3.5-45	22 Fr	3.5cm	45cm	7-10ml	1 Each

AVANOS* MIC* GASTRO-ENTERIC FEEDING TUBES/KITS

Kit Includes:

1- Gastro-Enteric Tube

1- .025in. Guidewire (for 16 Fr, 18 Fr Tubes)

1- .035in. Guidewire (for 20 Fr, 22 Fr, 24 Fr, 26 Fr, 28 Fr, 30 Fr Tubes)

Legacy Stock #	Size (Outer Dia.)	Jejunal Length	Balloon Volume	Packaging
0210-16LV	16 Fr	25.4cm	3-5ml Balloon	1 Each
0210-18LV	18 Fr	25.4cm	3-5ml Balloon	1 Each
0210-16	16 Fr	57.9cm	7-10ml Balloon	1 Each
0210-18	18 Fr	57.9cm	7-10ml Balloon	1 Each
0210-20	20 Fr	57.9cm	7-10ml Balloon	1 Each
0210-22	22 Fr	57.9cm	7-10ml Balloon	1 Each
0210-24	24 Fr	57.9cm	7-10ml Balloon	1 Each
0210-26	26 Fr	57.9cm	7-10ml Balloon	1 Each
0210-28	28 Fr	57.9cm	7-10ml Balloon	1 Each
0210-30	30 Fr	57.9cm	7-10ml Balloon	1 Each

+Jejunal length measured from the distal aspect of the retention balloon to the catheter tip.

ACCESSORIES

AVANOS* MIC-KEY* FEEDING TUBE ACCESSORIES

For use with MIC-KEY* G, J and GJ Feeding Tubes.
DEHP-Free Formulation [Di-(2-Ethylhexyl)phthalate]

Legacy Stock #	ENFit® Stock #	Length	Description	Packaging
0121-12	0141-12	12in.	MIC-KEY* Extension Set with SECUR-LOK* Right-Angle Connector and 2 Port "Y" and Clamp	5 Each
0121-24	0141-24	24in.	MIC-KEY* Extension Set with SECUR-LOK* Right-Angle Connector and 2 Port "Y" and Clamp	5 Each
0122-02	0142-02	2in.	MIC-KEY* Medication Set with SECUR-LOK* Right-Angle Connector and 2 Port "Y"	5 Each
0123-12	0143-12	12in.	MIC-KEY* Bolus Extension Set with Catheter Tip, SECUR-LOK* Straight Connector and Clamp	5 Each
0123-24	0143-24	24in.	MIC-KEY* Bolus Extension Set with Catheter Tip, SECUR-LOK* Straight Connector and Clamp	5 Each
0124-12	N/A	12in.	MIC-KEY* Bolus Extension Set with Catheter Tip, SECUR-LOK* Right-Angle Connector and Clamp	5 Each
0124-24	N/A	24in.	MIC-KEY* Bolus Extension Set with Catheter Tip, SECUR-LOK* Right-Angle Connector and Clamp	5 Each
98460	N/A		MIC-KEY* Stoma Measuring Device	10 Each

AVANOS* MIC* JEJUNOSTOMY FEEDING TUBE ACCESSORIES

For Use With the MIC* Jejunostomy Feeding Tube

Legacy Stock #	Description	Packaging
0301-50	Jejunostomy Repair Kit – Includes Silicone Med A Medical Adhesive, Stent, Syringe, Blunt Tip Needle, Suture, Silicone Overcollar	1 Each
0301-55	Jejunostomy Luer-Lok® Caps	5 Each
0301-60	Jejunostomy Feeding Port Adapter	5 Each

AVANOS* MIC* FEEDING TUBE ACCESSORIES

For Use with MIC* Gastrostomy, Bolus Gastrostomy, PEG, Jejunal, Gastro-Enteric, and GJ Feeding Tubes. DEHP-Free Formulation [Di-(2-Ethylhexyl)phthalate]

Legacy Stock #	Length	Description	Packaging
0105-06	6in.	Extension Tubing with Bolus and Stepped Connectors at Opposite Ends	1 Each
0105-12	12in.	Extension Tubing with Bolus and Stepped Connectors at Opposite Ends	1 Each

AVANOS* MIC* PERCUTANEOUS ENDOSCOPIC GASTROSTOMY FEEDING TUBE ACCESSORIES

For Use With MIC* Percutaneous Endoscopic Gastrostomy Feeding Tubes

Legacy Stock #	ENFit® Stock #	Description	Packaging
0135-14	8135-14	Universal Feeding Adapter – For use with 14 Fr MIC* PEG	1 Each
0135-20	8135-20	Universal Feeding Adapter – For use with 20 Fr MIC* PEG	1 Each
0135-24	8135-24	Universal Feeding Adapter – For use with 24 Fr MIC* PEG	1 Each
0136-14	8136-14	Bolus Feeding Adapter – For use with 14 Fr MIC* PEG	1 Each
0136-20	8136-20	Bolus Feeding Adapter – For use with 20 Fr MIC* PEG	1 Each
0136-24	8136-24	Bolus Feeding Adapter – For use with 24 Fr MIC* PEG	1 Each

AVANOS* TRANSITION CONNECTORS

ENFit® Stock #	Description	Packaging
7149-00	Distal End ENFit® Transition Connector	200 Qty/Case
F00071	Distal End ENFit® Transition Connector	30 Qty/Case
F00073	Oral Syringe Connector: Allows an Oral Syringe to be connected to an ENFit® feeding tube	30 Qty/Case
F00074	Transition Funnel Connector: Allows an ENFit® feeding tube to accept a Step ("Christmas Tree") Connector or a Catheter Tip Syringe	30 Qty/Case

NOTES:

AVANOS

avanosmedicaldevices.com



*Registered Trademark or Trademark of Avanos Medical, Inc., or its affiliates.
© 2018 AVNS. All rights reserved.
EC172000