



## OUR MISSION:

To become the most trusted name in electrosurgery.

## OUR FOCUS:

Provide innovative solutions that improve patient safety and outcome, are simple and easy to use, and offer superior value.

## OUR COMMITMENT:

Deliver unrivaled service and safe and effective medical devices that meet regulatory requirements and exceed customer expectations.

## OUR LEGACY:


To be recognized as a company of integrity that does right by our customers, employees, communities and the environment.

To order call: 1-800-747-6110 +1-801-576-9669









## Table of Contents

Corporate Overview	4
E-Z Clean® Non-Stick Electrosurgical Tips	7
Coated Blade Electrodes	9
Coated Needle Electrodes	10
Blunt Needle Electrodes	10
Coated Ball Electrodes	11
Coated Specialty Electrodes	11
Protective Nose Cone E-Z Clean Electrodes	13
MEGAfine® Needle Electrodes	15
ACE Blade™ Advanced Cutting System	17
Pencils	21
E-Z Pen™ Reusable Electrosurgical Pencils	23
Disposable Pencils	25
Patient Return Electrodes	27
Mega Soft® & Mega Soft Dual Cord	29
Mega 2000®	29
Mega Soft Cables	30
Disposable Patient Return Electrodes	31
Laparoscopic Electrodes	33
E-Z Clean® Disposable Lap Electrodes	35
Accessory	36
E-Z Clean® Disposable MEGAtips®	38
Reusable Indicator Shafts®	39
Accessory	39
Reusable Laparoscopic Electrodes	41
Accessory	41
All-In-One® Hand Control Device	43
All-In-One E-Z Clean® Laparoscopic Electrodes	43
Electrosurgical Generator	45
Mega Power™ Generator	47
Accessories	47
Smoke Evacuation	49
Mega Vac™ Smoke Evacuators	51
Smoke Evacuator Pencils and Shrouds	51
Accessories	52
Suction Coagulators	53
Lletz Loop Electrodes	57
Square Loop	59
Round Loops	59
Bipolar Forceps	61
Educational Products Summary	65

Megadyne Catalog-Rev e2b.indd 2-3 1/27/10 10:19 A

## **MEGADYNE**®

The Electrosurgical Authority™

The Megadyne family of quality electrosurgical products utilizes new and innovative technology to improve safety, increase efficiency, and bring cost savings.
















# Megadyne—Your Trusted Partner in Electrosurgery.



**MEGADYNE:** The electrosurgical products company was founded in 1985 with the launch of the E-Z Clean® line of non-stick, eschar resistant, PTFE coated electrodes and a mission of facilitating the research and development of high quality, cost effective medical products for surgeons and nurses.

After more than 25 years of innovation, Megadyne's research and development team continues to develop and bring to market innovative solutions that improve patient safety and outcome, are simple and easy to use, and offer superior value. The company's commitment to service, our customers, communities, employees and environment has established Megadyne as the Trusted Electrosurgical Authority.

Megadyne maintains a quality management system which conforms to the requirements of 21 CFR Part 820 (FDA Quality System Regulations), and is certified to ISO 13485, Canadian Medical Device Requirements, and Japanese MHLW Ordinance No. 169. In addition, devices we manufacture are certified to meet the requirements of the European Medical Device Directive 93/42/EEC and bear the CE mark.

**PRODUCTS:** Megadyne's products have been proven in thousands of operating rooms worldwide. From the Mega Power electrosurgical generator with proprietary ACE mode, to the patented Mega Soft patient return electrode, to the proprietary non-stick E-Z Clean pencils and electrodes, Megadyne provides innovative solutions that improve patient safety and outcome, are simple and easy to use, and offer superior value.

**CUSTOMERS:** Megadyne's innovative products have been used in millions of electrosurgical procedures worldwide. In addition to hospitals and surgical centers, our customers include healthcare group purchasing organizations and supply companies in over 75 countries throughout the world.



- E-Z Clean® Non-Stick Electrodes
- ACE Cutting System
- E-Z Pen™ Electrosurgical Pencils
- Disposable Electrosurgical Pencils
- MEGAfine® Needle Electrodes
- Mega Soft® Patient Return Electrodes
- The Mega Power<sup>™</sup>
   Electrosurgical Generator
- Mega Vac<sup>™</sup> Smoke Evacuation System
- Disposable Laparoscopic Electrodes
- Resposable Laparoscopic Electrodes
- Reusable Stainless Steel Laparoscopic Electrodes
- All-In-One® Hand Control
- LEEP/LLETZ Loop Electrodes
- Suction Coagulators
- Bipolar Forceps
- Electrosurgical Cables

# Why Megadyne is "The Trusted Electrosurgical Partner"

### PRODUCT INNOVATIONS

- Over 25 years experience developing electrosurgical devices.
- 1st to develop the non-stick tip. E-Z Clean® tips continue to be the most preferred non-stick tip in the market today. E-Z Clean® allows the surgeon to minimize any delays due to eschar build-up and to use the lowest possible power settings to achieve the desired effect, which minimizes collateral thermal necrosis.
- 1st to develop a safe reusable patient return electrode that maintains safe current densities. The Mega Soft line is rapidly replacing the 30+ year old "sticky" pad.
- 1st to develop a non-adhesive patient return electrode for pediatric patients between 0.3 and 50 pounds which protects their fragile skin and maintains safe current densities.
- 1st to develop a patient return electrode that incorporates pressure reduction.
- 1st to develop the Advanced Cutting Effect (ACE)
   System which provides the precision of a scalpel with
   the power and speed of electrosurgery. The ACE System
   saves time and reduces/eliminates the use of a scalpel
   in the sterile field.
- 1st to develop the reusable laparoscopic Indicator Shaft® which provides increased safety from damaged insulation and reduces cost.

- 1st to develop the reusable electrosurgical pencil with a reliable counting mechanism and end-of life indicator. The E-Z PEN™ provides a reliable, cost effective pencil with an E-Z Clean® tip at an "everyprocedure-use" price.
- A very intuitive and easy-to-use electrosurgical generator. Mega Power™ is designed as a fully functioning generator that benchmarked the best of the best and incorporated those core features into a simple, easy to use device. This is truly an innovation, and is offered at a competitive price.
- Continuing education programs available to customers at any time.
- "Preventing Surgical Fires", provides one nursing contact hour in a lecture program presented by our Surgical Specialists.
- "Laparoscopic Electrosurgical Safety", provides one nursing contact hour in a lecture type program presented by our Surgical Specialists.
- "It's so Easy Being Green", provides one nursing contact hour in a web-based format.
- "Basics of Electrosurgery", a lecture presented by our Surgical Specialists (No continuing education).
- Through innovative and superior technological advances in reusable products, Megadyne provides cost savings to the OR. Megadyne delivers the industry's best patient safety, clinical benefits, continuing education and cost effectiveness.

Megadyne's continued commitment to research and development provides clinicians with the safest, highest quality products available from the electrosurgical industry's complete solution provider. From the Mega Power electrosurgical generator to active electrodes to return electrodes, Megadyne provides you with the complete circuit of electrosurgical patient care.



**MEGADYNE E-Z**CLEAN E-Z Clean Non-Stick Electrosurgical Tips E-Z Clean® non-stick electrodes feature a proprietary polytetrafluoroethylene (PTFE) coating that reduces eschar buildup during surgical procedures, enabling surgeons to use lower power settings. Lower power settings mean less thermal necrosis to surrounding tissue and a further reduction in eschar build-up Cost-Effective Non-Stick Electrosurgery Tips For a Variety of Surgical Procedures

## E-Z CLEAN NON-STICK ELECTROSURGICAL TIPS

**E-Z**CLEAN°



E-Z CLEAN NON-STICK
ELECTROSURGICAL TIPS

## Trusted by Thousands of Surgeons Worldwide



- **New:** modified tips are insulated everywhere but the distal 3-5 mm of the electrode shaft with PTFE insulation that resists splitting in temperatures beyond 700° F, minimizing the likelihood of damage to surrounding tissue.
- E-Z Clean tips utilize a pure PTFE coating that resists flaking at high temperatures unlike competitive tips that often utilize cheap polyurethane and other fillers in their PTFE coatings.
- With more than two-dozen tip configurations and sizes, Megadyne offers blade, needle, and ball electrodes.
- Surgeons can get the performance of Megadyne's E-Z Clean non-stick tips for every procedure at costs comparable to stainless steel with our pencil and tip savings program. Contact us at 800-747-6110 for a free sample and to learn about this program.

IF IT ISN'T GREEN, IT ISN'T MEGADYNE E-Z CLEAN®

Description	Catalog Number	Qty Per Cas
Coated Blade Electrodes		
MEGADYNE		
2.5 inch (6.35 cm), Standard Blade Electrode	0012	12
MEGADYNE		
2.5 inch (6.35 cm), Standard Blade Electrode Modified	0012M	12
MEGADYNE DE LA FILLA DE LA FIL	2012	40
2.75 inch (7 cm), X-Long Standard Blade Electrode	0012A	12
MEGADYNE		
2.75 inch (7 cm), X-Long Standard Blade Electrode Modified	0012AM	12
<b>FZCLEAN</b>		
4 inch (10.2 cm), Extended Blade Electrode	0014A	12
MEGADYNE ZZCLEAN		
4 inch (10.2 cm), Extended Blade Electrode Modified	0014AM	12
MEGADYNE' EZICLEAN'		
6.5 inch (16.5 cm), Extended Blade Electrode	0014	12
MEGADYNE ZZCLEAN		
6.5 inch (16.5 cm), Extended Blade Electrode Modified	0014M	12

# WITE BY WITE

## E-Z CLEAN NON-STICK ELECTROSURGICAL TIPS

**MEGADYNE E-Z**CLEAN°

Description	Catalog Number	Qty Per Case
Coated Needle Electrodes		
MEGADYNE		
2.75 inch (7 cm), Needle Electrode	0013	12
MEGADYNE		
2.75 inch (7 cm), Needle Electrode Modified	0013M	12
MEGADYNE EZCLEAN		
4 inch (10.2 cm), Extended Needle Electrode	0016A	12
MEGADYNE' 154CLEAN		
4 inch (10.2 cm), Needle Electrode Modified	0016AM	12
MEGADYNE EZICLEAN		
6 inch (15.2 cm), Extended Needle Electrode	0016	12
MEGADYNE #ZCLEAN	•	
6 inch (15.2 cm), Extended Needle Electrode Modified	0016M	12

## Blunt Needle Electrodes

MEGADYNE		
2.75 inch (7 cm), Blunt Needle Electrode	0113	12
MEGADYNE		
2.75 inch (7 cm), Blunt Needle Electrode Modified	0113M	12
MEGADYNE E-ZCLEAN		
4 Inch (10.2 cm), Blunt Needle Electrode	0113A	12

Description	Catalog Number	Qty Per Cas
Coated Ball Electrodes		
MECADIONE:		
2 inch (5.1 cm), Ball Electrode	0015	12
MEGADYNE ZZCLEAN		
5 inch (12.7 cm), Ball Electrode	0009	12
Coated Specialty Electrodes		
MEGADYNE 15	ZCLEAN	_
Bayonet Needle Electrode	0028	12
MEGADYNE E	CLEAN	
Bayonet Needle Electrode Modified	0028M	12
MEGADYNE E	CLEAN	

## E-Z CLEAN NON-STICK ELECTROSURGICAL TIPS

MEGADYNE E-ZCLEAN®



PROTECTIVE NOSE
CONE ELECTRODES

# FOR HIGH-MOISTURE OR CONFINED-SPACE PROCEDURES

- New protective nose cone electrodes provide patient safety in procedures performed in moisture-rich areas and confined spaces by minimizing the risk of alternate site injuries.
- Protective nose cone reduces risk of inadvertent thermal burns caused by moisture-induced alternate current pathways
- Minimizes chances of moisture wicking into the hand piece
- Features new modified PTFE insulation that resists splitting in temperatures beyond 700° F
- Proprietary E-Z Clean® PTFE-coated tips save time by reducing eschar build-up

IF IT ISN'T GREEN, IT ISN'T MEGADYNE E-Z CLEAN®

Description	Catalog Number	Qty Per Case
Protective Nose Cone E-Z Clean® Electro	odes	
E-ZICLEAN		
2.5 inch (6.35 cm), Blade Electrode Modified with PTFE Insulation and Nose Co	ne 0012MD	24
DZ CLEAN		
2.75 inch (7 cm), Blade Electrode Modified with PTFE Insulation and Nose Cone	e 0012AMD	24
[E-ZISLEAN]		
4 inch (10.2 cm), Blade Electrode Modified with PTFE Insulation and Nose Cone	e 0014AMD	24
IE-ZIGIBANI		
6.5 inch (16.5 cm), Blade Electrode Modified with PTFE Insulation and Nose Col	ne 0014MD	24
E-ZICLEAN		

## E-Z CLEAN NON-STICK ELECTROSURGICAL TIPS

MEGADYNE MEGAFINE°




NEEDLE ELECTRODES

MEGAfine®
Needle Electrodes
Enable Precise
Electrosugery at Lower
Power Settings



- Very precise MEGAfine tip allows for lower power settings which help to control thermal necrosis and damage to surrounding tissues
- Less eschar buildup for less cleaning and a faster case
- Up to 60% cost savings compared to competitors
- MEGAfine extra-sharp needles are coated with the proprietary E-Z Clean PTFE coating

IF IT ISN'T GREEN, IT ISN'T MEGADYNE E-Z CLEAN®



## **MEGADYNE®**

## E-Z CLEAN NON-STICK ELECTROSURGICAL TIPS

NOTES:
OI ETE CA
COMPTENT
MEGADYNE
TOSURGICAL

**MEGADYNE ACE**BLADE<sup>™</sup> **ACE: Advanced Cutting Effect** System This single, innovative instrument—combined with the ACE Mode of the Mega Power Generator—effectively cuts, coagulates and dissects, reducing the need to pass scalpels and conventional electrosurgical blades back and forth. ACE Blade glides through tissue and safely performs skin incisions with wound healing equivalent to a scalpel. The Precision of a Scalpel with the Power and Speed of Electrosurgery

## •

## ACE ADVANCED CUTTING SYSTEM

MEGADYNE ACEBLADE\*\*



ACE ADVANCED
CUTTING SYSTEM

The Precision of a Scalpel with the Power and Speed of Electrosurgery



CE SKIN INCISION



## Time Savings

- The single, innovative instrument combined with the proprietary ACE Mode of the Mega Power Generator, can be safely used from skin to skin: incision, coagulation and dissection
- With its unique geometry, ACE Blade glides through tough tissue with ease and precision
- E-Z Clean® PTFE non-stick coating minimizes eschar buildup, saving valuable time



## **Improved Safety & Efficacy**

- Reduces or eliminates the need for a scalpel in the OR
- Cuts like a scalpel without having a razor's edge
- Produces wound healing equivalent to a scalpel



## **Superior Control**

- Smoother and straighter incisions without requiring tension
- Easier to create vertical skin edges for challenging incisions

Description	Catalog Number	Qty Per Case
ACE Blade Electrodes		
ACE BLADE		
ACE Blade 2.75" (7 cm)	ACE12A	12
ACE BLADE		
ACE Blade 4" (10.2 cm)	ACE14A	12
AGE BLADE		
ACE Blade 6.5" (16.5 cm)	ACE14	12
Pencils with ACE Blades*		
ALEGADYNE § OF S	)	
Rocker Switch Pencil with ACE Blade and Holster	ACE30H	50
ALEGADINE \$ 0 S	1	
Button Switch Pencil with ACE Blade and Holster	ACE35H	50


Mega Power™ Generator with ACE Mode 1000



\*Available in the U.S. Only


## **MEGADYNE**°

## ACE ADVANCED CUTTING SYSTEM



20 To order call: 1-800-747-6110

# MEGADYNE E-ZCLEAN®



## PENCILS



MEGADYNE E-ZPEN™



E-Z PEN REUSABLE
ELECTROSURGICAL PENCIL

# It's E-Z to Save Time and Reduce Operating Costs

Megadyne's new cost-per-procedure pricing offers surgical facilities reusable pencils and E-Z Clean blades at a low per-procedure cost. By partnering with Megadyne, hospitals can provide surgeons with the E-Z Clean nonstick tips they prefer for the same price or less than they are currently paying.

- The safety lockout mechanism prevents usage after 12 uses.
- 12-use lifespan promotes cost efficiencies
- Lightweight, ergonomic design with kink resistant cord for comfortable handling
- E-Z Clean PTFE coating lasts longer and requires less cleaning than other coatings
- Water resistant construction minimizes
   electrical shocks and inadvertent activations
- Push Button or Rocker Switch configurations

Description Catalog Number Qty Per Case

E-Z Pen Reusable Electrosurgical Pencil







Button Switch Pencil with 10 ft. (3.0 m) cord (with holster)
Button Switch Pencil with 15 ft. (4.5 m) cord (with holster)


3035H 10 3035HL 10

## •

## **PENCILS**



PENCILS



## DISPOSABLE PENCILS

Whether You Prefer a Rocker Switch, Or Like To Push Buttons, Megadyne Has an Electrosurgical Pencil For You.



- Ergonomic and lightweight design is comfortable and minimizes fatigue
- Water resistant construction minimizes electrical shocks and inadvertent activations
- Cords are flexible and memory-free for easier handling of pencil
- Packaged with an E-Z Clean or stainless steel blade

Description Catalog Number Qty Per Case

## Disposable Pencils

MEGADYNE	30 00 00	MEGADYNE
Rocker Switch Pencil with E-Z Clean Electrosurgical Tip	0030	50
Rocker Switch Pencil with E-Z Clean Electrosurgical Tip and Holster	0030H	50

MEGADYNE	TO OOGE	MEGADYNE
Button Switch Pencil with E-Z Clean Electrosurgical Tip	0035	50
Button Switch Pencil with E-Z Clean Electrosurgical Tip and Holster	0035H	50

E-Z	SO S	MEGADYNE
Rocker Switch Pencil with Modified E-Z Clean Electrosurgical Tip	0036	50
Rocker Switch Pencil with Modified E-Z Clean Electrosurgical Tip and Holster	0036H	50

	0046	MEGADYNE
Button Switch Pencil with Modified E-Z Clean Electrosurgical Tip	0037	50
Button Switch Pencil with Modified E-Z Clean Electrosurgical Tip and Holster	0037H	50

## Pencils with ACE Blades\*

SQVI8 EDVI	- S S S S S S S S S S S S S S S S S S S	MEGADYNE	
Rocker Switch Pencil with ACE Blade and Holster	ACE30H	50	



\*Available in the U.S. only

To order call: 1-800-747-6110 | 25



Megadyne Catalog-Rev e2b.indd 24-25

## **MEGADYNE®**

## PENCILS

Description	Catalog Number	Qty Per Case
Disposable Pencils		

		8 NECADI	RE
Rocker Switch Pencil with Stainless Steel Electrosurgical Tip	0038	50	
Rocker Switch Pencil with Stainless Steel Electrosurgical Tip and Holster	0038H	50	



Button Switch Pencil with Stainless Steel Electrosurgical Tip	0039	50
Button Switch Pencil with Stainless Steel Electrosurgical Tip and Holster	0039H	50



PENCIL HOLSTER

# MEGADYNE MEGADYNE

Patient Return Electrodes

The Mega Soft eliminates the small disposable sticky pad that can damage patients' skin.

If the patient is in contact with only a small portion of the pad, safety still prevails because the Mega Soft is a self-contained current limiting device. This innovative technology makes Mega Soft the 'Guardian of Electrosurgical Patient Safety'.


Safety and Cost Savings in Electrosurgery



## PATIENT RETURN ELECTRODES

MEGADYNE MEGASOFT®

35 Million+ Procedures. Zero Pad-Site Burns.




RETURN ELECTRODES

## IT'S OFFICIAL: **MEGA SOFT IS SAFER THAN REM\***







Excessive heating and pad-site burns using "sticky pad" CQM return electrode

- Safer than sticky-style monitoring pads. Large electrode size and capacitive "gatekeeper" design ensures patient safety.
- Easy set-up. No shaving or searching for an adequate patient pad site. "As simple as making the bed".
- Kinder to the patient. No adhesive-related skin irritations and all the benefits of a pressure reduction pad.
- Pediatric Mega Soft specially designed for patients between 0.8 and 50 lbs. (.36-23 Kg)
- Environment friendly. Eliminates the waste of disposable sticky pads.
- Saves money. The reusable Mega Soft significantly reduces cost per procedure—up to a 70% savings.
- Confidence. We provide a warranty and holdharmless for each pad
- Radiolucent and latex-free.

Description

**Catalog Number** 

**Qty Per Case** 

## Mega Soft Pediatric, Mega Soft & Mega Soft Dual Cord



Mega Soft Pediatric Patient Return Electrode 66 cm (26") L x 30.5 cm (12") W, 1.3 cm (1/2") thickness



Mega Soft Patient Return Electrode 117 CM (46") L x 51 cm (20") W, 1.25 cm (1/2") thickness



Mega Soft Dual Cord Patient Return Electrode 117 CM (46") L x 51 cm (20") W, 1.25 cm (1/2") thickness 0835

## Mega 2000®

Mega 2000 Reusable Patient Return Electrode	(36" L x 20" W-91 cm x 51 cm)	0800	1
Mega 2000 Sheath		0825	30

To order call: 1-800-747-6110 29

REM is a trademark of ValleyLab

## PATIENT RETURN ELECTRODES

Description	Catalog Number	Qty Per Case
Mega Soft & Mega Soft Dual Cord Cables		
Mega Soft Reusable DetachaCable™ Cord, Standard Connector, 2.4m (8.0 ft.)	M2K-01	1
Mega Soft Reusable DetachaCable™ Cord, Standard Connector, 4.4m (14.4 ft.)	M2K-02	1
Mega Soft Reusable DetachaCable™ Cord, Phone Plug Connector, 2.4m (8.0 ft.)	M2K-03	1
Mega Soft Reusable DetachaCable™ Cord, Phone Plug Connector, 4.4m (14.4 ft.)	M2K-04	1
Mega Soft Reusable DetachaCable™ Cord, Extended Phone Plug Connector, 4.4m (14.4 ft.)	M2K-05	1
Mega Soft Reusable DetachaCable™ Cord, Argon Beam Connector, 2.4m (8.0 ft.)	M2K-06	1
Replacement Pigtail Cable	M2K-07	1
DetachaCable™ Cord, Compatibility, Mega Soft 2.4m (8.0 ft.)	M2K-08	1
DetachaCable™ Cord, Compatibility, Mega Soft 4.4m (14.4 ft.)	M2K-09	1
Mega Soft Patch Kit	PKITO01	1

30 To order call: 1-800-747-6110

## PATIENT RETURN ELECTRODES



DISPOSABLE

Megadyne disposable patient return electrodes with disposable or reusable cables expand the options available to you.



- Compatible with most electrosurgical generators
- Low-profile flexible backing easily conforms to patient contours
- Enhanced adhesive border and gel overlap of foil edge eliminate foil exposure and risk of incidence
- Return electrode contact quality monitoring style pads minimize the potential for pad-site burns

To order call: 1-800-747-6110 31

Megadyne Catalog-Rev e2b.indd 30-31 1/27/10 10:22 AM

## PATIENT RETURN ELECTRODES

Description	Catalog Number	Qty Per Case
Pre-Attached Cords		
Adult, single plate, with 3m (10') pre-attached cord	0850C	50
Pediatric single plate, with 3m (10') pre-attached cord, 2.7-11.4 kg (6-25 lbs.)	0860C	50
Adult, dual plate, with 3m (10') pre-attached cord	0855C	50
Pediatric dual plate, with 3m (10') pre-attached cord, 2.7-11.4 kg (6-25 lbs.)	0865C	50
ReCordable™		
Adult, single plate, no cord	0850	50
Pediatric single plate, no cord, 2.7-11.4 kg (6-25 lbs.)	0860	50
Adult, dual plate, no cord	0855	50
Pediatric dual plate, no cord, 2.7-11.4 kg (6-25 lbs.)	0865	50
Reusable Cable, Single Style, 3m (10')	0870	1
Reusable Cable, Dual Style, 3m (10')	0875	1

32 To order call: 1-800-747-6110

# MEGADYNE E-ZCLEAN®

# Laparoscopic Electrodes

E-Z Clean electrosurgical instruments
provide excellent dissection and
coagulation for a variety of surgical
applications, rendering procedures cleaner,
safer and faster. They have a patented,
proprietary, non-stick surface which
resists eschar buildup and last an entire
case—saving both time and money.

E-Z Clean is the Laparoscopic Electrode of Choice for Thousands of Surgeons.





## LAPAROSCOPIC ELECTRODES

**E-Z**CLEAN°



Coated Laparoscopic
Electrodes Save
Surgeons Valuable Time
in the Operating Room.

- High-tech E-Z Clean Coating for Reduced Eschar Buildup
- Provide Excellent Dissection and Coagulation for a Variety of Surgical Applications
- Reintroduction During Surgery is Minimized
- Lower Power Settings for Safer Delivery of Electrosurgery
- Disposable Electrodes provide a new fully insulated shaft for each case, preventing stray current and unwanted burns

IF IT ISN'T GREEN, IT ISN'T MEGADYNE E-Z CLEAN®

Description Standard Cat # Split Stem Cat # Qty Per Case

Disposable Laparoscopic Electrodes Compatible With Standard Pencils

Compatible With Reusable Foot Control Cords

MEGADYNE ZZCLEAN 6

MEGADYNE ZICLEAN

33 cm. Spatula

0018

6

 MEGADYNE EZICLEAN

 33 cm. Curved
 0019
 0019S
 6

 45 cm. Curved
 0019L
 0019LS
 6

 MEGADYNE EZICIEAN

 33 cm. L-Hook
 0020
 0020S
 6

 45 cm. L-Hook
 0020L
 0020LS
 6

 MEGADYNE EZICLEAN

 33 cm. J-Hook
 0021
 0021S
 6

 45 cm. J-Hook
 0021L
 0021LS
 6

 MEGADYNE
 ₱₹/CLEAN

 33 cm. Needle
 0022
 0022S
 6

Continued...

To order call: 1-800-747-6110 35

Megadyne Catalog-Rev e2b.indd 34-35



## LAPAROSCOPIC ELECTRODES

Description	Standard Cat #	Split Stem Cat #	Qty Per Case
Disposable Laparoscopic Electrodes	Compatible With Standard Pencils	Compatible With Reusa Foot Control Cords	ble

	MEGADYNE ZCLEAN	
33 cm. Ball	0024	6

	MEGADYNE 152CLEAN		
33 cm. O-Wire	0026S	6	_

	MEGAD	YNE <b>3</b> 2CLEAN	
33 cm. L-Wire	0100	0100S	6
45 cm. L-Wire	0100L	0100LS	6

		MEGADYNE #ZCLEAN	
33 cm. J-Wire		0105S	6
Accessory	STANDARD STEM (FOR ELECTROSURGICAL PENCIL)	SPLIT STEM (FOR FO	OT CONTROL CORD USE)


0075

3 m (10-foot) 1

36 To order call: 1-800-747-6110

## LAPAROSCOPIC ELECTRODES





- Indicator Technology Minimizes Risk and Promotes Safety Megadyne's patented inner yellow "Indicator" insulation technology alerts
- clinicians when nicks or holes in the outer insulation occur indicating it is time to discard and replace the electrode.
- E-Z Clean Tips for Less Eschar Buildup. Megadyne's High-tech PTFE coating reduces eschar buildup and easily wipes clean with a sterile sponge making procedures cleaner, faster and safer.
- Reusable Shafts for Cost Savings E-Z Clean MEGAtips are cost-effective when combined with Reusable Indicator Shafts.

To order call: 1-800-747-6110 37

Reusable Foot Control Cord

## LAPAROSCOPIC ELECTRODES

MEGATIP\*

Description	Catalog Number	Qty Per Case
Disposable MEGAtips		
MEGADYNE		
L - Wire Disposable MEGAtip	0600	6
MEGADANE		
L - Wire Modified Disposable MEGAtip	0600M	6
C MEGADYNE		
J - Wire Disposable MEGAtip	0605	6
MEGANDYNE		
Curved Spatula Disposable MEGAtip	0618	6
· MEGADYNE		
Curved Blade Disposable MEGAtip	0619	6
MEGADYNE		
L - Hook Disposable MEGAtip	0620	6
MEGADYNE		
L - Hook Modified Disposable MEGAtip	0620M	6
MEGADYNE		
J - Hook Disposable MEGAtip	0621	6
MEGADYNE		
J - Hook Modified Disposable MEGAtip	0621M	6

Description	Catalog Number	Qty Per Case
Reusable Indicator Shafts		
32 cm Reusable Indicator Shaft	0690	2
32 cm Reusable Indicator Shaft with Split Stem	0690S	2
38 cm Reusable Indicator Shaft	0695	2
38 cm Reusable Indicator Shaft with Split Stem	0695S	2
Accessory	SPLIT STEM (FOR FOOT	CONTROL CORD USE)
Reusable Foot Control Cord	0075	3 m (10-foot)

Full product indemnification and warranty are available. See your Megadyne representative for details.

## LAPAROSCOPIC ELECTRODES

MEGADYNE INDICATORSHAFT\*



REUSABLE
LAPAROSCOPIC ELECTRODES

Megadyne—the leader in electrosurgical innovation—offers a new technology for safety-conscious surgeons.

Indicator Technology Minimizes
 Risk and Promotes Safety

Megadyne's patented inner yellow "Indicator" insulation technology alerts clinicians when nicks or holes in the outer insulation occur indicating it is time to discard and replace the electrode.

• A Collection of Tip Configurations

Our stainless steel laparoscopic electrodes are available in five tip configurations including Spatula, Curved, L-Hook, J-Hook, and L-Wire.

• Control and Surgical Precision

Indicator Shafts are sturdy and rigid, allowing direct tactile control of the tip, and mastery over tough fibrotic tissue and blunt dissection.

INDICATOR SHAFT

DescriptionStandard Cat #Split Stem Cat #Qty Per CaseReusable Laparoscopic ElectrodesCompatible With Standard PencilsCompatible With Reusable Foot Control Cords

Reusable Indicator Curved Spatula Laparoscopic Tip - 32 cm 3018 3018S 2
Reusable Indicator Curved Spatula Laparoscopic Tip - 38 cm 3018A 3018SA 2

Reusable Indicator Curved Laparoscopic Tip - 32 cm 3019 3019S 2
Reusable Indicator Curved Laparoscopic Tip - 38 cm 3019A 3019SA 2

Reusable Indicator L-Hook Laparoscopic Tip - 32 cm 3020 3020S 2
Reusable Indicator L-Hook Laparoscopic Tip - 38 cm 3020A 3020SA 2

Reusable Indicator J-Hook Laparoscopic Tip - 32 cm 3021 3021S 2
Reusable Indicator J-Hook Laparoscopic Tip - 38 cm 3021A 3021SA 2

Reusable Indicator L-Wire Laparoscopic Tip – 32 cm 3100 3100S 2
Reusable Indicator L-Wire Laparoscopic Tip – 38 cm 3100A 3100SA 2

STANDARD STEM

Accessory

CONTROL CORD USE)

Reusable Foot Control Cord

0075

3 m (10-foot)

Full product indemnification and warranty are available. See your Megadyne representative for details.

To order call: 1-800-747-6110 41

SPLIT STEM (FOR FOOT

## LAPAROSCOPIC ELECTRODES

MEGADYNE ALL-IN-ONE®



ALL-IN-ONE®
HAND CONTROL DEVICE

Megadyne's All-In-One Hand Control provides complete laparoscopic control at your fingertips: Irrigation, aspiration and electrosurgery.



- High-tech PTFE coating for reduced eschar build-up
- Convenient electrosurgery with suction and irrigation in one hand control device
- Lower power settings for safer delivery of electrosurgery
- Eschar wipes clean with a sterile sponge saving time in the O.R.
- Sliding electrode sheath provides pinpoint suction

Description	Catalog Number	Qty Per Cas
All-In-One Hand Control Device		
MASABANAS		
All-In-One Hand Piece	0055	5
All-In-One Coated Laparoscopic Ele	ectrodes	
L-Hook - 32cm working length	0060	5
J Wire-Hook - 38cm working length	0085	5
Curved - 32cm working length	0062	5
'L' Wire Hook - 32cm working length 'L' Wire Hook - 38cm working length	0063 0080	5 5
Standard - 32cm working length	0064	5
All-In-One Electrosurgical Tip		
	0066	

## **MEGADYNE**°

## LAPAROSCOPIC ELECTRODES

NOTES:	
	MPLETE CIRC
	MEGADYNE
	TOSURGICAL

44 To order call: 1-800-747-6110

MEGADYNE MEGAPOWER<sup>™</sup>
Mega Power<sup>™</sup>

Electrosurgical Generator

In addition to offering the most advanced electrosurgical tools in the industry, Megadyne also provides a state-of-the art Mega Power electrosurgical generator—designed for its simplicity and ease of use, as well as its powerful and elegant features.

Simple on the outside, advanced on the inside.



MEGADYNE

MEGA ROWER

## MEGA POWER ELECTROSURGICAL GENERATOR

MEGADYNE

**MEGA** POWER™

**Catalog Number** 



MEGA POWER
ELECTROSURGICAL GENERATOR

## The Power of Simplicity and Performance



- Mega Power's proprietary Advanced Cutting Effect (ACE) mode—in combination with the ACE Blade—is designed to provide a scalpel-like cutting effect for minimal thermal necrosis and reduced scarring.
- Large, easy to read displays that are highly visible, even at a distance.
- The Constant Control Technology™ feature monitors tissue impedance and adjusts power output to reduce tissue damage and drag.
   The result is a smooth, clean, accurate cutting effect at the lowest possible settings for maximum patient safety.
- The recall feature readily recalls previous power settings and modes on set-up.

Mega Power Generator



Mega Power Generator 1000

## Accessories



ega Cart 1300



art Shelf 1300S



Monopolar Footswitch (3 meter cable) 1400



Bipolar Footswitch (3 meter cable) 1450



Monopolar Adapter 96007

Mega Power's standard warranty is a full two years. An 8 year extended warranty is available for Megadyne partners. See your Megadyne representative for details.

## **MEGADYNE®**


## MEGA POWER ELECTROSURGICAL GENERATOR

**NOTES:** 

48 To order call: 1-800-747-6110

## **MEGADYNE MEGA** VAC

# Smoke **Evacuators**






FLOW SET

## MEGA VAC SMOKE EVACUATOR

MEGADYNE MEGA VAC



SMOKE EVACUATOR PENCILS

Protecting you and your patient from the hazards of electrosurgical smoke plume is easy, quiet and comfortable with the Mega Vac smoke evacuator system.



- Variable Levels of Suction—turning on and off when the surgeon activates and deactivates the electrosurgical pencil based upon the surgeon's discretion.
- Use-Time Monitoring—saving filter expense by letting the user know precisely when to change the filter.
- Creates effective suction with less noise during electrosurgical procedures.
- Ultra Vac Pencil and Attacha Vac Pencil Shroud allow Surgeons to easily capture smoke while comfortably performing surgery.

Description **Catalog Number Qty Per Case** Smoke Evacuators



MEGA Vac™ Smoke Evacuator	2100	1	
MEGA Vac Plus™ Smoke Evacuator with Laparoscopic Mode	2200	1	

## Smoke Evacuator Pencils and Shrouds



To order call: 1-800-747-6110 | 51

Megadyne Catalog-Rev e2b.indd 50-51

## MEGA VAC SMOKE EVACUATOR

Description	Catalog Number	Qty Per Case
Connectors	-	
Connector, Filter/Tubing	2140	1
Connector, Surgimedics, some Buffalo Filter 22mm Male to 10mm Male	2145	1
Connector, Stackhouse, Niche, some Buffalo Filter 22mm Female to 10mm Male	2155	1
Accessories		
Laparoscopic Tubing	2190	25
Speculum Tubing	2195	25
ULPA Filter with Fluid Trap	2210	9
THE PROPERTY OF THE PROPERTY O		
Charcoal Output Filter	2220	1
RF Sensor	2250	1
Manual Foot Switch	2160	1

MEGADYNE SUCTION COAGULATORS

Suction Coagulators

These disposable suction coagulators provide precise coagulation with controlled fluid evacuation. The ergonomic handle combined with a flexible, kink-resistant cannula maximizes suction while allowing the surgeon to customize the instrument to the specific procedure at hand.

Uncompromising Standards of Safety



# ®

## SUCTION COAGULATORS

## MEGADYNE SUCTION COAGULATORS



# Designed With the ENT Surgeon In Mind, This Line of Ergonomic Suction Coagulators Provides:

- Variety of sizes to meet all patient parameters
- Comes in hand or foot controlled version for surgeon's preference
- Fully insulated for patient safety
- Highest quality at a competitive price
- Integrated "Shock protector" keeps surgeon's finger safe
- Disposable, to reduce risk of CJD transmission

Description	French Sizes	Catalog Number	Qty Per Case
Suction Coagulators			
MEGADYNE			
Foot Control, 6" (15.2 cm) cannula, 10ft. (3.0 m) cord	10	0040-10	10
Foot Control, 6" (15.2 cm) cannula, 10ft. (3.0 m) cord	10	0040-25	25
AUSGARDYNIS			
Hand Control, 6" (15.2 cm) cannula, 10ft. (3.0 m) cord	10	0041-25	25
Hand Control, 6" (15.2 cm) cannula, 10ft. (3.0 m) cord	8	0042-25	25
Hand Control, 6" (15.2 cm) cannula, 10ft. (3.0 m) cord	12	0043-25	25



## **MEGADYNE®**

**NOTES:** 

## SUCTION COAGULATORS

OMPLETE CIAC
MEGADYNE
TO STATE OF THE ST
USURGICA

56 To order call: 1-800-747-6110

# MEGADYNE LLETZLOOP

# Lletz Loop Electrodes


These disposable Lletz Loop Electrodes are
easy to use and will provide a precise specimen
for examination. Each Lletz Loop Electrode
configuration delivers excellent cutting
characteristics for precise tissue excision.

Uncompromising Standards of Safety



## LLETZ LOOP ELECTRODES

MEGADYNE LLETZLOOP



LLETZ LOOP


# Let Us Put You Into Our Loop For Cost Savings...High-Quality At a Competitive Price

- Color-coded for easy size and shape recognition
- High-quality at a competitive price
- Disposable for user convenience
- Fine tungsten wire loop for excellent 'cutting' characteristics
- Stable "Y" construction to minimize loop flex
- Fits most standard electrosurgical pencils

Description	Catalog Number	Qty Per Case
Square Loop		
MEGADYNE		
10 mm Wide x 7 mm Deep with 11.8 cm Shaft Orange	0410	12
Round Loops		
MEGADYNE		
10 mm Wide x 5 mm Deep with 11.8 cm Shaft Purple	0420	12
MEGADYNE		
15 mm Wide x 6 mm Deep with 11.8 cm Shaft Gray	0430	12
MEGADYNE		
20 mm Wide x 10 mm Deep with 11.8 cm Shaft Red	0440	12
MEGADYNE		
10 mm Wide x 10 mm Deep with 11.8 cm Shaft Yellow	0450	12
MEGADYNE		
15 mm Wide x 12 mm Deep with 11.8 cm Shaft Green	0460	12
		continued.

## LLETZ LOOP ELECTRODES

Description	Catalog Number	Qty Per Case
Round Loops continued		
MEGADYNE		
20 mm Wide x 12 mm Deep with 11.8 cm Shaft White	0470	12
MEGADYNE		
20 mm Wide x 15 mm Deep with 11.8 cm Shaft Blue	0480	12
MEGADYNE SOURCE CONTROL OF THE STATE OF THE	0.400	10
20 mm Wide x 8 mm Deep with 11.8 cm Shaft Tan	0490	12



## BIPOLAR FORCEPS

MEGADYNE BIPOLARFORCEPS



BIPOLAR FORCEPS
AND ACCESSORIES

# Megadyne's bipolar forceps provide the latest technology available including the following features:

- Highly polished tips reduce tissue adherence
- Light closing force causes less fatigue over time
- Well balanced, fitting comfortably in the surgeon's hand
- Nylon coating resists abrasion and provides superb insulation
- Precision tip alignment aids fine control
- Reusable, durable, and designed for multiple uses

Description	Catalog Number	Qty Per Case
Bipolar Forceps		
		_
		_
Scoville-Greenwood Forcep, 1.5mm tip, 7 in. (17.8 cm)	4010	1
ACCULATION AND ACCULATION	_	
Jeweler Forcep, 0.4mm tip, 4 in. (10.2 cm)	4020	1
Cushing forcep, 1.5mm tip, 6 in. (15.2 cm)	4030	1
Adson forcep, 1.0mm tip, 4 in. (10.2 cm)	4040	1

continued...



## BIPOLAR FORCEPS

Description	Catalog Number	Qty Per Case
Bipolar Forceps		
· · · · · · · · · · · · · · · · · · ·		
Semkin forcep, 0.5 mm tip, 5 in. (12.7 cm)	4050	1
Gerald Forcep, 1.0 mm tip, 7 in. (17.8 cm)	4060	
anni 📗		
Bipolar Cable, disposable, flying lead, 12 ft. (3.7 m)	4000	10
Fixed Lead Cable	4005	10

## **MEGADYNE**

# Educational Products Summary



## EDUCATIONAL PRODUCTS SUMMARY

## **Continuing Education**

Megadyne is committed to customer training and education. Please contact your local Surgical Specialist to take part in the following continuing education and formal in-service programs, and to receive your copy of the educational booklet, Principles of Electrosurgery.



Preventing Surgical Fires—PowerPoint Presentation

• 1 contact hour program, lecture



Laparoscopic Electrosurgical Safety—PowerPoint Presentation

• 1 contact hour program, lecture



It's so Easy Being Green (Your Guide to a Greener O.R.)


• 1 contact hour, web-based

www.megadyne.com/education



Principles of Electrosurgery booklet

Handout, no continuing education




The Basics of Electrosurgery—PowerPoint Presentation

• Lecture, no continuing education







1506 South State Street | Draper, Utah 84020 | Phone (801) 576-9669 | Toll-Free 800-747-6110 | www.megadyne.com

To order call: 1-800-747-6110 +1-801-576-9669

©Copyright 2010 Megadyne Medical Products, Inc.

MMP 3000119-01 Rev. E 2010-01-2