









Product Family
Patient Eyeshields - Disposable

We offer a range of patented disposable eye patches for use with Lasers, IPL (Intense Pulsed Light), LEDs (Light Emitting Diodes), Low-level Laser Therapy treatments and Microdermabrasion.



Product Numbers & Ordering Information

Product Numbers	Details
 31-7300	<p>Derm-Aid, 50 pr/box Provides a physical barrier for use during common dermatological procedures such as microdermabrasion, chemical peels and facials. CE marked.</p>
 31-7303	<p>Derm-Aid, Sample Pack, 5 pair Provides a physical barrier for use during common dermatological procedures such as microdermabrasion, chemical peels and facials. CE marked.</p>
 31-8300	<p>Laser-Aid, 24 pr/box Provides protection from all lasers operating between 190-10,600nm. These shields also disperse the heat typically emitted from infrared lasers. Protection Level: OD 7 @ 190-11,000nm, EN 207: 315-1400 D L6 + IR L7, >1400-11000 D L4 + IR L3. CE marked.</p>
 31-8303	<p>Laser-Aid, Sample Pack, 5 pair Provides protection from all lasers operating between 190-10,600nm. These shields also disperse the heat typically emitted from infrared lasers. Protection Level: OD 7 @ 190-11,000nm, EN 207: 315-1400 D L6 + IR L7, >1400-11000 D L4 + IR L3. CE marked.</p>
 31-8400	<p>LED Aid, 50 pr/box Provides protection from LED devices of all wavelengths from 315-1400nm and low-powered (<1 watt) lasers. Blocks out the bright light from red, yellow and blue LED's and IPL systems. Protection Level: OD 4 @ 315-1400nm, EN 207: >315-1400 DIR L4. CE marked.</p>
 31-8403	<p>LED Aid, Sample Pack, 5 pair Provides protection from LED devices of all wavelengths from 315-1400nm and low-powered (<1 watt) lasers. Blocks out the bright light from red, yellow and blue LED's and IPL systems. Protection Level: OD 4 @ 315-1400nm, EN 207: >315-1400 DIR L4. CE marked.</p>

Patient Eyeshields - Disposable



31-9300

IPL-Aid, 50 pr/box

Provides protection from the wavelengths emitted from IPL devices hence blocking out the bright light commonly associated with IPL systems. CE marked.



31-9303

IPL-Aid, Sample Pack, 5 pair

Provides protection from the wavelengths emitted from IPL devices hence blocking out the bright light commonly associated with IPL systems. CE marked.



31-7500

Anesthesia-Aid, 50 pr/box

Protect the eye during anesthesia by keeping the eye moist, prevent corneal abrasions and keep particles out of the eyes. CE marked.



31-7503

Anesthesia-Aid, Sample Pack, 5 pair

Protect the eye during anesthesia by keeping the eye moist, prevent corneal abrasions and keep particles out of the eyes. CE marked.

Hazardous Materials

None

Overview

Key Features

These Sperian eye shields enable full access to temples, eyebrows and nose-bridge and are designed to protect most orbital sizes. Thanks to these innovative products, cleaning and disinfecting patient or client eyewear is a thing of the past! They are simple and quick to use, latex free and hypoallergenic. They remain securely in place during the procedure. The convenient tabs make them easy and painless to remove or reposition.

- Allows access to nose bridge, brow and temples
- Easy to apply and comfortable to wear
- Tab for easy removal or repositioning
- Easy-peel adhesive
- Latex-free and hypoallergenic

Laser-AID Protection Level:

- OD 7 @ 190-11000nm
- EN 207: 315-1400 D L6 + IR L7, >1400-11000 D L4 + IR L3

LED-AID Protection Level:

- OD 4 @ 315-1400nm
- EN 207: >315-1400 DIR L4

Derm-AID, IPL-AID:

- CE marked.

Recommended Industries/Use

- Medical

Hazards

- Laser/IPL Light

Regulations

- ANSI Z136.1 - Meets ANSI Z136.1 standard
- EN207 - Meets CE EN207. Laser-AID Protection Level: OD 7 @ 190-11000nm, EN 207: 315-1400 D L6 + IR L7, >1400-11000 D L4 + IR L3; LED-AID Protection Level: OD 4 @ 315-1400nm, EN 207: >315-1400 DIR L4, Derm-AID, IPL-AID: CE marked.

-
-
-

Historical Brand

Glendale

Specifications

Standards

- Meets Standards: ANSI Z136.1, CE EN207.Laser-AID, Protection Level: OD 7 @ 190-11,000nm EN 207: 315-1400 D L6 + IR L7, >1400-11000 D L4 + IR L3.LED-AID, Protection Level: OD 4 @ 315-1400nm EN 207: >315-1400 DIR L4Derm_AID, IPL-AID: CE marked.

Literature & Documents

Data Sheet Laser Eyesheids_Eng

http://www.honeywellsafety.com/Supplementary/Documents_and_Downloads/Eye_and_Face_Protection/Patient_Eyeshields/19293/1033.aspx

Data Sheet Laser Eyesheids_Eng

Training & Proper Use

Please refer to the Instruction Sheet provided with the product.

Additional Information

laserinfo@sperianprotection.com or call 1-800-500-4739 (US) / +49-6131-9066499 (Int'l)

Technical Support

laserinfo@sperianprotection.com or call 1-800-500-4739 (US) / +49-6131-9066499 (Int'l)

RELATED PRODUCTS