



NEOMED[®]

NeoDrape[®] by NeoMed, Inc.

NEOMED IS PLEASED TO PRESENT THE NEODRAPE PRODUCT LINE

The NeoDrape product line offers effective thermoregulation and barrier protection for patients in Labor and Delivery and NICU's during sterile procedures. Within the NeoDrape product offerings are two separate items:

- The NeoDrape (Item # ND) is a sterile wrap that is used to swaddle the newborn patient during delivery through stabilization to maintain proper thermoregulation while minimizing moisture loss.
- The NeoDrape Fenestrated (Item # ND-FS) complies with recommendations set forth in the NANN PICC Guidelines for Practice, 2nd edition Pg 25, step 17: "Place a sterile drape underneath and above the insertion area. Cover as much of the infant as can be safely done while ensuring the ability of adequate observation."

Below are the common features and benefits of the NeoDrape and NeoDrape Fenestrated (Item # ND and ND-FS):

- Complies to the 2005 NRP guidelines for "Temperature Control"
- Meets the recommendations set forth by the Golden Hour
- Does not impede radiant heat common with paper drapes
- Allows for clear visualization of the patient during sterile procedures
- Inhibits dangerous static electricity, a must in high oxygen environments
- Manufactured medical grade polyethylene
- Cost effective
- Sterile packaging

NeoDrape (Item # ND)

1




2



3



NeoDrape Fenestrated (Item # ND-FS)



Hospital Product List

Features and Benefits of NeoDrape Item # ND:

- Used in the delivery room
- Allows for coverage from the back of the head through the top of the head held firm by a cap to minimize rapid heat loss from the head as demonstrated by the Golden Hour
- Sterile packaging allows wrap to be placed in sterile field for immediate use during delivery
- Measures 24 X 24 inches

Features and Benefits of NeoDrape Fenestrated Item # ND-FS:

- Complies with new NANN PICC guidelines for Practice, 2nd edition: Pg 25, step 17
- Used during sterile procedures
- Fenestration allows for easy placement and removal over surgical site
- Measures 32 X 50 inches creating a large sterile field from side rail to side rail
- Has adhesive strips on both sides to allow for easy securement to drop side bassinet

Product Code	Product Description	Quantity
NeoDrape		
ND	NeoDrape Occlusive Drape/Wrap	50/case
ND-FS	NeoDrape Fenestrated and Perforated Drape/Wrap	40/case
Catheterization Tray		
CT	Catheterization Tray - No Catheter	5/Case
CT-NS	Catheterization Tray - No Catheter, No Saline	5/Case

NeoDrape Fenestrated is packaged inside our Catheterization Trays for UAC/UVC insertion

5-0811